

# PHILIPS

IPL – Hair removal device

Lumea Advanced

For use on body

15 minutes to treat lower legs

>250,000 lamp flashes

Extra long cord

SC1993/70



## Gentle, long-term hair reduction

Philips Lumea Advanced works wonders to prevent the reappearance of hair on body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

### Long lasting hair removal with IPL technology

- Enjoy smooth skin every day
- Adapted for safe and effective use in your own home
- Effortlessly effective
- Suitable for a wide variety of hair and skin types

### Safe, fast and effective on body

- For use on the body
- Large treatment window for fast application
- Gentle treatment

### Maximum freedom, minimum hassle

- No replacement parts, no hidden costs
- Five adjustable light energy settings


# Highlights

## IPL technology




Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

## Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. As a consequence, the amount of hair your body grows gradually decreases. Repeating the treatment leaves your skin beautifully hair free and touchably smooth.

## Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four fortnightly treatments resulting in smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth. – To enlarge this picture, please click on the image in the image gallery on the top of this page

## Suitable hair and skin types



Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blond or white/grey hair. This product is also NOT suitable for dark skin.

## Application areas



Use with confidence to remove hair on body areas below the neckline: legs, armpits, bikini area, arms and stomach.

## Fast treatment on legs

Big treatment window for fast application on large areas like legs.

## Gentle treatment


Gentle to your skin.

## No replacement parts needed



Philips Lumea is a complete solution that comes ready to work right out of the box. More importantly, it does not require any replacement lamps or gels.

## Five light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

# Specifications

**Application areas**

**Body areas:** Arms, Legs, Belly, Bikini, Armpits

**Safety and adjustable settings**

**Integrated UV filter:** Protects skin from UV light

**5 light energy settings:** Adjustable to your skin type

**Integrated safety system:** Prevent unintentional flashing

**Application mode**

**Corded/cordless use:** Corded

**Stamp and Flash:** For treatment on curvy areas

**Power**

**Voltage:** 100–240 V

**Technical specifications**

**High performance lamp:** Delivers <gt/> 250,000 flashes

**Attachments**

**Body attachment (4 cm<sup>2</sup>):** For use below the neckline

**Service**

**Warranty:** 2 years global warranty

**Items included**

**Instructions for use:** User manual

**Storage:** Pouch

**Adapter:** 24V / 1500mA

**Accessory:** Cleaning cloth

**Application time**

**Lower leg:** 8 min

**Armpit:** 1 min

**Bikini line:** 1 min

