# PHILIPS

IPL hair removal  
system

Lumea

For use on body and face

15 minutes to treat lower legs

Lifetime >100.000 light  
pulses

Extra long cord

SC1996/60

## The ideal way to prevent hair regrowth for smooth skin at home

Philips Lumea IPL works wonders to prevent the reappearance of hair on body and face. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

### Easy to use

- Slide and Flash mode for a more convenient application
- Adjustable to your skin type
- No replacement parts, no hidden costs

### Long lasting smoothness

- Enjoy smooth skin every day
- Effortlessly effective

### Hair removal with IPL technology

- Adapted for safe and effective use in your own home

### Safe, fast and effective on the body and face

- Five adjustable light energy settings
- Integrated skin tone sensor for extra safety
- Precision attachment for extra safety on face
- Body attachment: fast treatment (15 min. for 2 lower legs)

# Highlights

## Enjoy smooth skin every day


Philips Lumea applies gentle pulses of light to the hair root. This results in the hair shedding naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth everyday.

## Effortlessly effective

Our clinical studies have shown significant hair reduction in as little as four bi-weekly treatments resulting with smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

## 5 light energy settings

Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

## Application areas

Use with confidence to remove hair on all body areas. The big body attachment (4cm<sup>2</sup>) enables fast treatment on large areas like legs. The precision attachment with additional light filter provides extra protection for the use on upper lip, chin and cheeks.

### Big body attachment, 4cm<sup>2</sup>

Big body attachment enables fast treatment on large areas like legs.

### Slide and Flash mode

Main application mode enables sliding motion for more convenient and continuous application.

### No replacement parts needed

Philips Lumea is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels.

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Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for dark skin.

## IPL technology

Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

## Skin Tone Sensor

For extra safety the integrated Skin Tone Sensor measures the treated skin complexion at the beginning of each session and occasionally during the session. If it detects a skin tone that is too dark for Lumea Essential it will automatically stop emitting pulses.

# Specifications

## Power

**Voltage:** 100-240 V

## Service

**Warranty:** 2 years global warranty

## Technical specifications

**High performance lamp:** Delivers <gt/> 100,000 flashes

## Application areas

**Body areas:** Arms, Legs, Belly, Bikini, Armpits

**Face areas:** Chin, Sideburns, Upper lip

## Attachments

**Precision attachment (2cm2):** For use below the eyes

**Body attachment (4cm2):** For use below the neckline

## Items included

**Instructions for use:** Quick Start Guide, User manual

**Storage:** Pouch

**Adapter:** 24V / 1500mA

## Safety and adjustable settings

**Integrated UV filter:** Protects skin from UV light

**5 light energy settings:** Adjustable to your skin type

**Integrated safety system:** Prevent unintentional flashing

**Skin Tone sensor:** Detects your skin tone

## Application time

**Lower leg:** 8 min

**Slide and Flash:** For convenient application

**Armpit:** 1 min

**Bikini line:** 1 min

**Face areas:** 1 min

## Application mode

**Corded / cordless use:** Corded