

PHILIPS

IPL – Hair removal
device

Lumea Prestige

For use on body

15 minutes to treat lower legs

>250,000 lamp flashes

Cordless design

SC2005/30



Feel continuously hair free*

Philips Lumea Prestige is the most effective cordless IPL solution and works wonders to prevent the reappearance of hair on body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Long lasting hair removal with IPL technology

- Enjoy smooth skin every day
- Effortlessly effective
- Adapted for safe and effective use in your own home
- Suitable for a wide variety of hair and skin types

Gentle and effective, even on sensitive skin

- Safe and effective even on sensitive areas
- Large treatment window for fast application
- Gentle treatment

Maximum freedom, minimum hassle

- Cordless convenience
- Five adjustable light energy settings
- No replacement parts, no hidden costs
- Slide and Flash mode for a more convenient application
- Lumea Prestige comes in a luxurious pouch


Highlights

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. As a consequence, the amount of hair your body grows gradually decreases. Repeating the treatment leaves your skin beautifully hair free and touchably smooth.

Effortlessly effective

Clinical results after treatment: **example 1** (leg)



Clinical results after treatment: **example 2** (bikini line)



Our clinical studies have shown significant hair reduction in as little as four fortnightly treatments resulting in smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual

hair regrowth. - To enlarge this picture, please click on the image in the image gallery on the top of this page

Suitable hair and skin types




Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for very dark skin.

Application areas



Use with confidence to remove hair on body areas below the neckline: legs, underarms, bikini area, arms and belly.

Fast treatment on legs




Big treatment window for fast application on large areas like legs.

Gentle treatment



Gentle to your skin.

Cordless convenience



Lumea's unique cordless design is perfect for maximum reach and ease of use.

No replacement parts needed



Philips Lumea is a complete solution that comes ready to work right out of the box. More importantly, it does not require any replacement lamps or gels.

Slide and Flash mode



Main application mode enables sliding motion for a more convenient and continuous application.

Specifications

Application areas

Body areas: Arms, Legs, Belly, Bikini, Armpits

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light

5 light energy settings: Adjustable to your skin type

Integrated safety system: Prevent unintentional flashing

Application mode

Slide and Flash: For convenient application

Corded/cordless use: Cordless

Stamp and Flash: For treatment on curvy areas

Power

Charging time: 100 min

Voltage: 100-240 V

Running time: Minimum 370 flashes

Technical specifications

High performance lamp: Delivers <gt/> 250,000 flashes

Attachments

Body attachment (4 cm2): For use below the neckline

Service

Warranty: 2 years global warranty

Items included

Instructions for use: User manual

Charger: Multi voltage charger

Storage: Luxurious pouch

Accessory: Cleaning cloth

Application time

Lower leg: 8 min

Armpit: 1 min

Bikini line: 1 min

