

PHILIPS

IPL hair removal
system

Lumea

For use on body and face

15 minutes to treat lower legs

Lifetime >100.000 light
pulses

Extra long cord

SC1996/60



The ideal way to prevent hair regrowth for smooth skin at home

Philips Lumea IPL works wonders to prevent the reappearance of hair on body and face. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Easy to use

- Slide and Flash mode for a more convenient application
- Adjustable to your skin type
- No replacement parts, no hidden costs

Long lasting smoothness

- Enjoy smooth skin every day
- Effortlessly effective

Hair removal with IPL technology

- Adapted for safe and effective use in your own home

Safe, fast and effective on the body and face

- Five adjustable light energy settings
- Integrated skin tone sensor for extra safety
- Precision attachment for extra safety on face
- Body attachment: fast treatment (15 min. for 2 lower legs)

Highlights

Enjoy smooth skin every day




Philips Lumea applies gentle pulses of light to the hair root. This results in the hair shedding naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth everyday.

Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four bi-weekly treatments resulting with smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

5 light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

Application areas

Use with confidence to remove hair on all body areas. The big body attachment (4cm²) enables fast treatment on large areas like legs. The precision attachment with additional light filter provides extra protection for the use on upper lip, chin and cheeks.

Big body attachment, 4cm²



Big body attachment enables fast treatment on large areas like legs.

Slide and Flash mode



Main application mode enables sliding motion for more convenient and continuous application.

No replacement parts needed



Philips Lumea is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels.

Suitable hair and skin types

Geeignete Haar- und Hauttypen / Types de peau et de poils adaptés		Die integrierte Hauttonsensoren prüft die Belichtung in den dunklen Hauttönen ab und passt die Energie der Lichtpulse automatisch an Ihre Haut an. Bitte beachten Sie die Hinweise zur Anwendung.				
Hautfarbe / Couleur de la peau	Haarfarbe / Couleur des poils	1	2	3	4	5
1	1	✓	✓	✓	✓	✓
1	2	✓	✓	✓	✓	✓
1	3	✓	✓	✓	✓	✓
1	4	✓	✓	✓	✓	✓
1	5	✓	✓	✓	✓	✓
2	1	✓	✓	✓	✓	✓
2	2	✓	✓	✓	✓	✓
2	3	✓	✓	✓	✓	✓
2	4	✓	✓	✓	✓	✓
2	5	✓	✓	✓	✓	✓
3	1	✓	✓	✓	✓	✓
3	2	✓	✓	✓	✓	✓
3	3	✓	✓	✓	✓	✓
3	4	✓	✓	✓	✓	✓
3	5	✓	✓	✓	✓	✓
4	1	✓	✓	✓	✓	✓
4	2	✓	✓	✓	✓	✓
4	3	✓	✓	✓	✓	✓
4	4	✓	✓	✓	✓	✓
4	5	✓	✓	✓	✓	✓
5	1	✓	✓	✓	✓	✓
5	2	✓	✓	✓	✓	✓
5	3	✓	✓	✓	✓	✓
5	4	✓	✓	✓	✓	✓
5	5	✓	✓	✓	✓	✓

Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for dark skin.

IPL technology

Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

Skin Tone Sensor

For extra safety the integrated Skin Tone Sensor measures the treated skin complexion at the beginning of each session and occasionally during the session. If it detects a skin tone that is too dark for Lumea Essential it will automatically stop emitting pulses.

Specifications

Power

Voltage: 100-240 V

Service

Warranty: 2 years global warranty

Technical specifications

High performance lamp: Delivers <gt/> 100,000 flashes

Application areas

Body areas: Arms, Legs, Belly, Bikini, Armpits

Face areas: Chin, Sideburns, Upper lip

Attachments

Precision attachment (2cm²): For use below the eyes

Body attachment (4cm²): For use below the neckline

Items included

Instructions for use: Quick Start Guide, User manual

Storage: Pouch

Adapter: 24V / 1500mA

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light

5 light energy settings: Adjustable to your skin type

Integrated safety system: Prevent unintentional flashing

Skin Tone sensor: Detects your skin tone

Application time

Lower leg: 8 min

Slide and Flash: For convenient application

Armpit: 1 min

Bikini line: 1 min

Face areas: 1 min

Application mode

Corded / cordless use: Corded

