



Whiter, healthier teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. HealthyWhite is proven to remove everyday stains in just two weeks, with the regular use of the Clean and White Mode.

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- · Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- · Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle

Offers you a customised cleaning experience

- Easy-start programme to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2 minute mode for outstanding clean
- Sensitive Mode: Gentle teeth and gum cleaning



HX6781/09

Specifications

Technical specifications

· Battery: Rechargeable

Design and finishing

· Colour: Frost White

Items included

- · Charger: 1 multi-voltage charger
- Charger cover: 1 charger cover with brush head holder.
- Brush head colour rings: 3 colour rings for easy sharing

Ease of use

- Brush head system: Easy snap-on brush heads for optimal hygiene
- Brushing time: 42 2-minute brushings OR 3 weeks
- Battery indicator: Illuminated icon blinks when it's time to recharge
- · Handle: Slim ergonomic design and rubber grip for

easy handling.

Cleaning performance

- Speed: Up to 31,000 brush strokes per minute
- Brush heads: 1 ProResults standard brush head for an All-Around Clean
- Performance: Removes more plaque than Sonicare Essence/Sonicare CleanCare
- Health benefits: Helps improve gum health
- Whitening benefits: Can whiten teeth up to two shades in two weeks
- Timer: Smartimer and quadpacer
- Modes: 3 Modes: Clean, Sensitive and Clean and White

Service

· Warranty: 2-year limited warranty

Power

· Voltage: Multi-voltage charger

Highlights

Patented sonic technology

Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck

Uniquely angled brush head neck makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Naturally whiter teeth

Dynamic Cleaning Action and direct broader contact with each tooth have been proven to remove everyday stains and help keep teeth naturally whiter.

Safe and gentle

Your Sonicare is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Smartimer

Two-minute timer helps ensure Dental Professional recommended brushing time

Quadpacer

30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Easy-start programme

Gently increases power over the first 14 uses to ease into the Sonicare experience

Clean and White Mode

2 minutes of Clean mode with an additional 30 seconds of White mode to focus on your visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks, with the regular use of the Clean and White

Clean Mode

Clean Mode: Standard 2 minute mode for outstanding clean

Gentle teeth and gum cleaning

Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums



Issue date 2012-05-04

Version: 1.0.1

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

12 NC: 0000 000 00000

www.philips.com