FREQUENTLY ASKED QUESTIONS

My gums bleed when I use AirFloss. Is that normal?

Bleeding gums may be a sign of infection and a little bleeding can be normal when starting a new oral care routine, as you may be cleaning in areas not previously reached. However, if bleeding is excessive or does not stop within a few days of regular use of AirFloss, consult a dental professional.

How can I tell AirFloss is working?

AirFloss is designed to remove plaque from between teeth and help improve gum health. Plaque is not visible, so these signs of improved gum health show that AirFloss is working: gums may bleed less over time, gum tissue may look pink instead of red and the parts of your gums between teeth will look pointy, not round with inflammation. Ask your dental professional about the health of your gums at your next check-up.

Will AirFloss remove everything from between my teeth?

AirFloss is designed to remove plaque from between your teeth as part of your daily oral care routine, and it should remove most food particles. However, it is not designed to remove large or sticky food particles.







Register your product and get support at

www.philips.com/welcome



USING YOUR SONICARE AIRFLOSS



STEP 1: CHARGE

- Before first use, charge AirFloss for 24 hours.
- · The indicator light flashes green until it is fully charged. A yellow light means the device is extremely low on charge.
- · When fully charged, AirFloss should last for two weeks of use (14 uses).

ACTIVATION BUTTON

WATER RESERVOIR

STEP 3: PRIME

- Press the power button to turn on AirFloss.
- · Prime the device by pressing the activation button six times or until it sprays a fine mist.



POWER

NOTE: AirFloss is meant for use on the front/outside surface of the teeth only.







STEP 2: PREPARE

- · Attach the nozzle, and fill the reservoir with water or mouthwash.
- The reservoir holds enough liquid for 1 to 2 uses.



STEP 4: ACTIVATE

- Proper placement is key in producing results with AirFloss.
- · Place the guidance tip between the teeth at the gum line as shown in the image above. Apply gentle pressure to direct spray between teeth, close lips around nozzle to avoid unwanted spray and press the activation button.
- · Glide the guidance tip along the gum line until you feel it settle between the next two teeth, and press the button again. Repeat for all spaces between teeth, both upper and lower.