

Philips Sonicare ProResults
Compact sonic toothbrush
head

1-pack

ProResults

HX6021

Focused Clean

The smaller brush head size helps you easily manoeuvre the brush to focus on hard-to-reach areas.

Thorough clean feeling in your mouth

- Angled, radial-trimmed bristles expand the area covered
- Contoured bristles fit the natural shape of your teeth
- Brush head design maximises sonic motion

Easy to manoeuvre to focus on hard-to-reach areas

- Small brush head size helps you focus in tight places

Designed to optimise performance

- Reminder bristles ensure you are getting an effective clean

A brush head that fits multiple handles


- Snap-on brush head

PHILIPS
sonicare

sense **and** simplicity

Highlights

Snap-on brush head



Easy brush head placement and handle cleaning

Angled radial-trimmed bristles




The curved surface area created by angled, radial-trimmed bristles helps maintain contact with the surfaces of the teeth and gums as the brush head pivots, covering a larger surface area without increasing the size of the brush head.

Contoured bristles

Bristles are trimmed in a contoured pattern to fit the natural shape of teeth. The peaks and valleys along the length of the brush head are ergonomically designed to follow teeth topography and fit between teeth so you're less likely to miss places when you brush.

Maximise cleaning performance



Thoughtful design and precision engineering ensure that when this brush head is combined with the sonic motion and flexible brushing modes, you get a thorough clean.

Reminder bristles

Reminder bristles let you know when to replace the brush head. After three months of normal use bristles exhibit fatigue, and brush heads are less effective. Replace your brush head every 3 months.

Small brush head

The contoured bristle design on a smaller surface cleans as thoroughly as other Sonicare brush heads while allowing you to reach problem areas your dentist may have pointed out. This design helps you reach and focus on tight places where plaque and stains can hide.

Specifications

Ease of use

- Suitable for these models: FlexCare+, FlexCare, HealthyWhite, HydroClean, EasyClean

Technical specifications

- Operating time: For optimal results, a new brush head is recommended every three months.

Weight and dimensions

- Dimensions brush head packaging: 18.3 h x 6.3 w x 4.8 d cm
- Weight brush head packaging: 0.026 kg

Items Included

- Hygienic travel cap: 1 pcs

