



Philips Sonicare Essence
Rechargeable sonic
toothbrush

HX5581

Better plaque removal

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

Provides a superior clean

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

Proven to improve oral health

- Clinically proven safe and gentle

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Offers you a customised cleaning experience

- Easy-start programme to acclimate to the Sonicare experience

PHILIPS
sonicare

sense and simplicity

Specifications

Technical specifications

- Battery: Rechargeable

Design and finishing

- Colour: White

Items included

- Charger: 1 charger

Ease of use

- Brush head system: Easy screw-on brush heads
- Brushing time: 28 2-minute brushings OR 2 weeks
- Battery indicator: Light in power button blinks to indicate when to recharge

Cleaning performance

- Speed: Up to 31,000 brush strokes per minute
- Brush heads: 1 e Series standard brush head
- Performance: Better plaque removal than a manual toothbrush
- Health benefits: Helps improve gum health
- Whitening benefits: Helps remove tea, coffee and other stains
- Timer: Smartimer and quadpacer

Service

- Guaranteed: 2-year limited warranty

Power

- Voltage: 220-240 V, 50-60 Hz

Highlights

Patented sonic technology

Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck

Uniquely angled brush head neck makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Safe and gentle

Your Sonicare is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Smartimer

Two-minute timer automatically tracks brushing time and turns the toothbrush off after two full minutes of brushing is complete.

Quadpacer

30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Easy-start programme

Gently increases power over the first 14 uses to ease into the Sonicare experience



Issue date 2011-09-23

Version: 1.0.1

12 NC: 0000 000 00000

© 2011 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com