



Philips Sonicare FlexCare
Rechargeable sonic
toothbrush with UV brush
head Sanitizer

3 modes

2 brush heads 1 soft travel case

FlexCare



HX6932/55

Advanced Cleaning

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- UV brush head Sanitizer - Cleaner brush, better oral health

Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology
- Improves gum health in only two weeks

Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Three different modes for a better cleaning experience



PHILIPS
sonicare

sense and simplicity

Highlights

Patented sonic technology



Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck

The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

UV brush head Sanitizer

All-in-one sanitizing, charging and storage. The UV brush head Sanitizer helps keep your brush head free of bacteria. It kills up to 99% of bacteria and viruses.*

Naturally whiter teeth



The Dynamic Cleaning Action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Safe and gentle

Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Improves gum health

This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the

gumline for improved gum health in just two weeks. Giving a superior cleaning in between teeth and removes significantly more plaque to a manual toothbrush

Smartimer



2-minute timer on this toothbrush helps ensure dental professional recommended brushing time

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Easy-start program

Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Three different cleaning modes

The three modes include a clean mode for an outstanding clean in 2 minutes. A sensitive mode for gentle yet thorough cleaning of the gums and a massage mode to stimulate the gum.

Specifications

Technical specifications

- Battery: Rechargeable

Design and finishing

- Color: White with sky blue grip and chrome front panel

Items included

- Charger: 1 multi-voltage charger with UV brush head sanitizer and 1 multi-voltage travel charger
- Brush head color rings: 3 color rings for easy sharing
- Sanitizer: 1 UV brush head sanitizer

Ease of use

- Brush head system: Easy snap-on brush heads for optimal hygiene
- Brushing time: 42 2-minute brushings OR 3 weeks
- Battery indicator: Illuminated icon indicates battery life
- Handle: Slim ergonomic design and rubber grip for easy handling.

Cleaning performance

- Speed: Up to 31,000 brush strokes per minute
- Brush heads: 1 ProResults standard brush head for an All-Around Clean and 1 ProResults compact brush head for a Precision Clean
- Performance: Removes up to 83% of plaque in hard-to-reach areas
- Health benefits: Helps improve gum health
- Whitening benefits: Helps remove tea, coffee and other stains
- Timer: Smartimer and quadpacer
- Modes: 3 Modes: Clean, Sensitive, Massage and 2 Routines: Max Care and GoCare

Service

- Warranty: 2-year limited warranty

Power

- Voltage: Multi-voltage charger

