Philips Sonicare AirFloss

EN: User manual

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Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Use AirFloss with your favorite mouthwash as part of your oral care routine. In addition to brushing, AirFloss helps reduce the plaque in between teeth in order to improve your gum health.

Important

Read this user manual carefully before you use the appliance and save it for future reference.

Danger

- Keep the charger away from water. Do not place or store it over or near water contained in a bathtub, washbasin, sink etc. Do not immerse the charger in water or any other liquid. After cleaning, make sure the charger is completely dry before you connect it to the mains.

Warning

- The mains cord cannot be replaced. If the mains cord is damaged, discard the charger.
- Always have the charger replaced with one of the original type in order to avoid a hazard.

- If the appliance is damaged in any way (AirFloss nozzle, handle and/or charger), stop using it.
- This appliance contains no serviceable parts. If the appliance is damaged, contact the Consumer Care Centre in your country (see chapter 'Warranty and support').
- Do not use the charger outdoors or near heated surfaces.
- This appliance can be used by children aged from 8 years and above and by persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and if they understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children unless they are older than 8 and supervised.

Caution

- Consult your dentist before you use this appliance if you have had recent oral or gum surgery, suffer from an ongoing dental condition, had extensive dental work in the previous 2 months, or have suspected issues with your dental work (filings, crowns, etc.).
- Consult your dentist if excessive bleeding occurs after using this appliance or if bleeding continues to occur after 4 weeks of use.
- If you have medical concerns, consult your doctor before you use this appliance.
- This Philips appliance complies with the safety standards for electromagnetic devices. If you have a pacemaker or other implanted device, contact your physician or the appliance manufacturer of the implanted appliance prior to use.
- Do not use attachments other than the ones recommended by the manufacturer.
- Do not exert excessive pressure on the nozzle tip of the nozzle.

- To avoid physical injury, do not fill the reservoir with hot water. Do not put foreign objects in the reservoir.
- Only use the AirFloss for its intended purpose as described in this user manual.
- To avoid physical injury, do not direct the spray under the tongue or into the ear, nose, eye or other sensitive area.
- Do not clean the AirFloss nozzle, the handle, the charger, the water reservoir and its cover in the dishwasher or a microwave.
- To avoid an ingestion hazard, do not put isopropyl alcohol or other cleaning fluids in the reservoir.
- The AirFloss is a personal care appliance and is not intended for use on multiple patients in a dental practice or institution.

Electromagnetic fields (EMF)

This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

Your Philips AirFloss (Fig. 1)

- 1 Nozzle tip
- 2 AirFloss nozzle
- 3 Activation button
- 4 Reservoir and cover

10 English

- 5 Handle
- 6 Charge indicator
- 7 Power on/off button
- 8 Removable nozzle holder (select models)
- 9 Charger

Note: The contents of the box may vary based on the model purchased.

Charging

Charge the AirFloss for 24 hours before first use.

- 1 Plug the charger into an electrical outlet. (Fig. 2)
- 2 Place the handle on the charger. (Fig. 3)

The charge indicator flashes green to indicate that the appliance is charging. (Fig. 4)

Note: When it is fully charged, the AirFloss has an operating time of up to two weeks (or 14 uses).

Note: When the charge indicators flashes yellow, the battery is low and needs to be recharged (fewer than 3 uses left).

Preparing for use

- 1 Align the nozzle and press the nozzle down onto the handle ('click') (Fig. 5)

Note: We advise you to place the tip facing the front of the handle as the illustration shows, but the nozzle can be placed in the handle either way without affecting the function. You can try both options and use what works best for you in your oral care routine.

Note: Replace the nozzle every 6 months for optimal results. Also replace the nozzle if it becomes loose or no longer 'clicks' onto the handle.

Note: The AirFloss only supports AirFloss nozzles. Do not attempt to use non-AirFloss nozzles on the AirFloss handle.

- 2 Open the cover of the reservoir ('click') (Fig. 6).
- 3 Fill the reservoir with mouthwash or water. (Fig. 7)

Tip: For the best result and a fresher experience, we advise you to use mouthwash with your AirFloss.

Note: Do not use mouthwash that contains isopropyl myristate (e.g. Dentyl Active or high concentrations of essential oils, as it may damage the appliance. This includes all mouthwash sold in glass and ceramic packaging.

- 4 Close the cover of the reservoir by pressing ('click') (Fig. 8).
- 5 Press the power on/off button to switch on the appliance. (Fig. 10)
The green charge indicator lights up.
- 6 Press and hold the activation button until spray comes out of the tip of the nozzle. (Fig. 9)

Using the appliance

- 1 Make sure the AirFloss is switched on. (Fig. 10)
- 2 Place the tip between the teeth and apply gentle pressure. Slightly close your lips over the nozzle when it is in place, to avoid splashing. (Fig. 11)

Tip: Use your finger to press the activation button.

Note: The AirFloss is intended to be used on the front/outside surface of the teeth only.

- 3 Press the activation button to deliver a single burst of air and microdroplets of mouthwash or water between the teeth.
Alternatively, hold the activation button to deliver continuous bursts of mouthwash or water approximately every 1 second. (Fig. 12)
- 4 Slide the guidance tip along the gum line until you feel it settle between the next two teeth.
- 5 Continue this procedure for all spaces including behind your back teeth. (Fig. 13)
- 6 To switch off the AirFloss, press the on/off button again.

Automatic shut-off function

The automatic shut-off function automatically switches off the AirFloss if it has not been used for 4 minutes.

Cleaning

If the AirFloss becomes clogged, or to clean the inside of the AirFloss, fill the reservoir with warm water and press the activation button until the reservoir is empty.

Do not clean the nozzle, handle or charger in the dishwasher.

Do not use cleaning agents to clean the appliance. The formulation or strength of some agents may damage the appliance.

Handle and nozzle

- 1 Remove the AirFloss nozzle from the handle. Rinse the nozzle after each use to eliminate the residue left in it. (Fig. 14)
- 2 Open the reservoir and rinse it under a faucet to eliminate residue. (Fig. 15)

Tip: You can also use a cotton swab to wipe the reservoir and remove the remaining residue.

Note: Do not immerse the handle in water.

Note: Failure to remove residue from the nozzle or reservoir may result in an unhygienic condition.

- 3 Wipe the entire surface of the handle with a damp cloth. (Fig. 16)

Charger

- 1 Unplug the charger before you clean it.
- 2 Wipe the surface of the charger with a damp cloth. (Fig. 17)

Storage

If you are not going to use the AirFloss for an extended period of time, follow the steps below.

- 1 Open the cover of the reservoir and empty it.
- 2 Press the activation button until no more spray comes out of the nozzle.
- 3 Unplug the charger.
- 4 Clean the AirFloss nozzle, handle and charger. Refer to chapter 'Cleaning'.
- 5 Store the AirFloss in a cool, dry place away from direct sunlight.

Ordering accessories

To buy accessories or spare parts, visit

www.shop.philips.com/service or go to your Philips dealer. You can also contact the Philips Consumer Care Center in your country (see the international warranty leaflet for contact details).

Recycling

- This symbol means that this product shall not be disposed of with normal household waste (2012/19/EU) (Fig. 26).
- This symbol means that this product contains a built-in rechargeable battery which shall not be disposed of with normal household waste (Fig. 27) (2006/66/EC). We strongly advise you to take your product to an official collection point or a Philips service center to have a professional remove the rechargeable battery.
- Follow your country's rules for the separate collection of electrical and electronic products and rechargeable batteries. Correct disposal helps prevent negative consequences for the environment and human health.

Removing the rechargeable battery

Please note that this process is not reversible.

- 1 To deplete the battery, repeatedly press the activation button until the AirFloss no longer produces any bursts of air. (Fig. 18)
- 2 Remove the AirFloss nozzle from the handle. (Fig. 19)
- 3 Firmly grab the handle with one hand at the top and one hand on the bottom. (Fig. 20)
- 4 Snap the handle into two parts. (Fig. 21)

Note: This step requires a good amount of physical force. The top and bottom half of the handle should be separated once you snap the handle.

- 5 Remove the bottom half of the handle. (Fig. 22)
- 6 Cut all 6 wires on the internal component with scissors. (Fig. 23)
- 7 Completely pull apart the internal component from the top part of the handle. (Fig. 24)

14 English

- 8** The battery is located inside the top half of the handle. Remove the battery and dispose of it properly (not in household waste). (Fig. 25)

Warranty and support

If you need information or support, please visit

www.philips.com/support or read the international warranty leaflet.

Warranty restrictions

The terms of the international warranty do not cover the following:

- AirFloss nozzle
- Damage caused by use of unauthorized replacement parts.
- Damage caused by misuse, abuse, neglect, alterations or unauthorized repair.
- Normal wear and tear, including chips, scratches, abrasions, discoloration or fading.

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©2017 Koninklijke Philips N.V. (KPNV). All rights reserved. Philips and the Philips shield are trademarks of KPNV. Sonicare and the Sonicare logo are trademarks of Philips Oral Healthcare, LLC and/or KPNV.