

# PHILIPS

## sonicare

Standard sonic  
toothbrush heads

ProResults

4-pack

HX6014/35



## All-around Clean

### Sonicare electric toothbrush head

Replacement brush head for Sonicare electric toothbrush, with larger surface area and contoured bristles that thoroughly cleans and massages teeth and gums.

#### **Designed to optimise performance**

- Reminder bristles ensure you are getting an effective clean

#### **Thorough clean feeling in your mouth**

- Brush head design maximises sonic motion
- Contoured bristles fit the natural shape of your teeth
- Angled, radial-trimmed bristles expand the area covered

#### **A brush head that fits multiple handles**

- Snap-on brush head

## Highlights

### Snap-on brush head



Easy brush head placement and handle cleaning

### Reminder bristles



Reminder bristles let you know when to replace the brush head. After three months of normal use bristles exhibit fatigue, and brush heads are less effective. Replace your brush head every 3 months.

### Maximise cleaning performance



Thoughtful design and precision engineering ensure that when these Philips Sonicare replacement brushes are combined with the sonic motion and flexible brushing modes, you get a thorough clean.

### Contoured bristles



These Philips Sonicare electric toothbrush bristles are trimmed in a contoured pattern to fit the natural shape of teeth. The peaks and valleys along the length of the brush head are ergonomically designed to follow teeth topography and fit between teeth so you're less likely to miss places when you brush.

### Angled radial-trimmed bristles



The curved surface area created by angled, radial-trimmed bristles helps maintain contact with the surfaces of the teeth and gums as the brush head pivots, covering a larger surface area without increasing the size of the brush head.

## Specifications

### Weight and dimensions

Dimensions brush head packaging: 21.6 (H) x 8.4 (W) x 2.2 (D) cm

Weight brush head packaging: 0.07 kg

### Technical specifications

**Operating time:** For optimal results, a new brush head is recommended every three months.

**Replacement:** Reminder Bristles fade when a replacement is needed.

### Features

ProResults brush head

### Items Included

Hygienic travel cap: 4 pcs

### Additional features

Contoured brush head

### Ease of use

**Suitable for these models:** FlexCare+, FlexCare, HealthyWhite, HydroClean, EasyClean, DiamondClean

