

# PHILIPS

## sonicare

Brush Head -  
InterCare

InterCare

4-pack

Mini

HX9014/01



## Reaches deeper between teeth

A better clean between teeth for healthy gums

With advanced tufting technology and extra-long bristles, InterCare helps to remove up to 6x more plaque between teeth than a manual toothbrush and improve gum health in just two weeks.

### **Designed to optimize performance**

- Reminder bristles let you know when to use a new InterCare

### **Advanced InterCare brush head technology**

- Anchor-free tufting allows for a denser bristle pack
- Designed for interdental cleaning, gum care and polishing

### **Fits multiple handles**

- Works with any Philips Sonicare click-on toothbrush

### **Better plaque removal and gum care**

- Clinically proven to improve gum health in just two weeks
- Removes up to 6x more plaque between teeth vs manual brushes

# Highlights

## Quickly improves gum health

Clinical studies show our new interdental brush head helps reduce inflammation and bleeding in gums. This in turn, helps to prevent gingivitis and gum recession, improving your overall gum health in just two weeks.

## Anchor-free tufting

InterCare features advanced anchor-free tufting for 34% more bristles than traditional staple-based tufting.

## Innovative bristle design


InterCare has an innovative bristle design for targeted oral health care: interdental cleaning, gum care and polishing. The extra-long light green bristles sweep away plaque caught deep between teeth, the dark green bristles remove plaque trapped along the gum line and the white bristles polish away surface stains for naturally whiter teeth.

## Easy click-on brush head




The InterCare brush head clicks on and off your brush handle, making it easy to change. It fits a wide range of Sonicare toothbrush handles: DiamondClean, FlexCare Platinum, FlexCare+, FlexCare, HealthyWhite and EasyClean.

## Reaches deeper between teeth



Featuring extra-long bristles to reach deep between teeth and along the gum line, our new interdental brush head helps remove up to 6x more plaque between teeth than a manual toothbrush.

## Reminder bristles



Blue reminder bristles fade to white when it's time to change your brush head. Cleaning efficacy diminishes after several months of normal use due to bristle wear on a visible and microscopic level so you should replace your brush head every three months.

# Specifications

## Ease of use

**Suitable for these models:** FlexCare+, FlexCare, HealthyWhite, EasyClean, DiamondClean, FlexCare Platinum

**Brush head system:** Easily snaps on/off the handle

## Operating time

**When to change your brush head:** Replace every 3 months

