

PHILIPS

sonicare

AirFloss Pro/Ultra -
Trial

AirFloss Pro/Ultra

1 nozzle

2 oz. BreathRx mouth rinse



HX8341/11



A triple burst of better gingival health

Efficacy meets compliance for infrequent flossers

Philips Sonicare AirFloss Pro/Ultra**** is clinically proven as effective as floss in improving gingival health.* In fact, up to 97% of inconsistent flossers who tried AirFloss Pro/Ultra showed improved gum health.*

Improves patients' oral health

- Up to 99.9% plaque removal***
- Clinically proven as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

The easiest way to effectively clean between teeth

- An easy way to improve cleaning between teeth
- Takes just 60 seconds: Point, Press, Clean!
- Fresh and clean with mouthwash

Highlights

Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.***

Micro-droplet technology




Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Cleans in just 60 seconds



With the AirFloss Pro/Ultra, patients can clean their entire mouth in less than 60 seconds a day. Patients simply select their desired burst frequency (single, double or triple) and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh and clean with mouthwash




Patients can use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities between



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between patients' teeth.

High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your patient's liking with single, double, or triple bursts for each press of the activation button.

Specifications

Items included

AirFloss Pro/Ultra handle: 1
AirFloss Pro/Ultra nozzle: 1
Antibacterial mouth rinse: 2oz
Charger: 1

Design and finishing

Color: Gray

Cleaning performance

Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction

Performance: Removes up to 99.9% plaque**

Health benefits: Cleans entire mouth in 60 sec, Improves gum health in 2 weeks

Customizable settings: Auto-burst or manual burst, Single, double, triple burst

For best results: Change nozzle every 6 months

Ease of use

Handle: Slim ergonomic design

Battery indicator: Handle LEDs show charge level

Battery Life: 1-4 weeks*****

Nozzle system: Easily snaps on and off

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Charge time: 24 hours to full charge

Ingredients

Antibacterial mouth rinse: 0.075% CPC, Zytex

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty



* When used in conjunction with a manual toothbrush and antimicrobial rinse in patients with mild to moderate gingivitis. AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. For more information, please visit www.philips.com/airfloss/faq.
* *** AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.
* ***** depending on burst setting used