

PHILIPS

sonicare

Rechargeable
Sonicare AirFloss

AirFloss

2 Handles 4 Nozzles

2 travel cases



HX8154



An easier way to floss

Gently and effectively improves oral health

If flossing isn't a regular part of your oral care routine, you're not reaching the bacteria that can build up between teeth causing infection. Sonicare AirFloss helps you to get a deep clean everyday by taking the hassle out of flossing.

Achieve a 30-second clean

- Simple one-button use; as easy as a mouse-click

Simple: point & press

- Microburst technology cleans at the touch of a button
- Charging after two weeks
- Slim, angled nozzle for easy to even the hard-to-reach areas
- Guidance tip for easy placement
- Easy-to-fill reservoir
- Slim, easy to hold handle
- Fits standard Sonicare chargers

Proven improvement of oral health

- Up to 99.9% plaque removal


Highlights

Microburst technology



Microburst technology applies a quick burst of pressurized air and micro-water droplets to clean deep between teeth where a toothbrush can't reach

Slim, angled nozzle



Slim, angled nozzle for easy to even the hard-to-reach areas

Slim, easy to hold handle



Easy to manoeuvre around the mouth.

Charging after two weeks



The AirFloss lasts two weeks between charging.

Guidance tip



Simply glide the guidance tip along the gum-line until you feel it settle between the teeth. Provides easy positioning even on those hard-to-reach back teeth.

Fits Sonicare chargers



You only need one charger to charge your Sonicare products.

Simple one-button use




Simple one-button use for easy cleaning in 60 seconds.

Easy-to-fill reservoir



Easy-to-fill reservoir uses only a tiny amount of water for less mess. Fill with mouthwash for enhanced oral health benefits.

Up to 99.9% Plaque Removal



AirFloss uses microdroplets of air and water or mouthwash to remove up to 99% more plaque in-between teeth*

Specifications

Ease of use

Handle: Slim ergonomic design

Battery indicator: Light in power button blinks to indicate when to recharge

Technical specifications

Battery: Rechargeable

Charging time: 24 hour(s)

Operating time (full to empty): 2 weeks or 14 full-mouth uses

Operating time: Replace nozzle every 6 months

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty

