

HX8211/06



### Don't floss? Then AirFloss.

### Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

#### Achieve a 30-second clean

- · Cleans your entire mouth in just 30 seconds
- · Point. Press. Clean!

#### Helps build a healthy habit

 ${}^{\textstyle \cdot}{}$  An easy way to improve cleaning between teeth

#### Improves oral health

- · Removes more plaque between teeth
- ${}^{\raisebox{.5ex}{$\raisebox{.5ex}{}}}}}}}}}}}}}}}$
- Helps prevent cavities between teeth

#### Easy to use

- · Air and micro-droplet technology
- Comfortable ergonomic handle
- · Guidance tip ensures correct placement



## Highlights

#### Air and micro-droplet technology



Air and micro-droplet technology

#### An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers\* continued to use AirFloss four or more days per week.

#### **30 Second Clean**



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

#### Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks than manual brushing alone.

Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

#### **Guidance tip**



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth

#### Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

#### Removes more plaque between teeth



Removes more plaque between teeth

#### **Ergonomic handle**

Designed to be easy to maneuver in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.



#### Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

Ease of use

Handle: Slim ergonomic design

Battery indicator: Flashing indicates low

battery

Battery Life: Up to 3 wks when fully charged

Nozzle system: Easily snaps on and off

**Technical specifications** 

Battery: Rechargeable

Operating time: Replace nozzle every 6

months

Battery type: Lithium ION

Charge time: 24 hours to full charge

**Power** 

Voltage: Multi-voltage charger

**Service** 

Warranty: 2-year limited warranty

**Items included** 

AirFloss Handle: 1 AirFloss Nozzle: 1

Charger: 1 multi-voltage charger

Nozzle holder: 1

**Cleaning performance** 

Performance: Removes up to 5x more plaque in-between teeth than manual brushing alone Health benefits: Improves gum health in only two weeks than manual brushing alone

Auto-burst feature: Cleans entire mouth in 30

sec

Nozzles: Best results change every 6mos

**Design specifications** 

Color(s): White with green accents

© 2019 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2019-07-23 Version: 2.0.1

EAN: 08 71010 36607 29

www.philips.com

