

PHILIPS

sonicare

AirFloss Pro/Ultra -
Interdental cleaner

w/ 1 nozzle

2 oz. BreathRx mouth rinse



HX8331/12



Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra**** is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.**

Improves oral health

- Up to 99.9% plaque removal***
- Clinically proven as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology


- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

The easiest way to effectively clean between teeth

- An easy way to start a healthy habit
- Takes just 60 seconds: Point, Press, Clean!
- Total daily confidence of a fresh and clean mouth


Highlights

Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.***

Micro-droplet technology




Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh & clean with mouthwash




Fill the reservoir on the handle with either mouthwash or water, then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

Specifications

Items included

AirFloss Pro/Ultra handle: 1

AirFloss Pro/Ultra nozzle: 1

Charger: 1

Mouth rinse: 2oz BreathRx

Design and finishing

Color: White with pink accents

Cleaning performance

Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction

Performance: Removes up to 99.9% plaque**

Health benefits: Cleans entire mouth in 60 sec, Improves gum health in 2 weeks

Customizable settings: Auto-burst or manual burst, Single, double, triple burst

For best results: Change nozzle every 6 months

Ease of use

Handle: Slim ergonomic design

Battery indicator: Handle LEDs show charge level

Battery Life: 1-4 weeks*****

Nozzle system: Easily snaps on and off

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Charge time: 24 hours to full charge

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty



* or your money back

* * when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.

* * * In a lab study, actual in-mouth results may vary.

* * * * AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.

* * * * * depending on burst setting used