

PHILIPS

sonicare

AirFloss Pro/Ultra -
Interdental nozzles

3 nozzles

HX8033/26



Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro nozzle is the easiest way to effectively clean between teeth.

Improves oral health

- Clinically proven as effective as floss for gum health**
- Up to 99.9% plaque removal***
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

The easiest way to effectively clean between teeth

- An easy way to start a healthy habit


Highlights

Better gum health




Philips Sonicare AirFloss Pro/Ultra nozzle is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Triple burst customization




AirFloss Pro nozzle's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.***

High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro nozzle helps prevent cavities from forming in the spaces between your teeth.

Specifications

Items included

AirFloss Ultra nozzle: 3

Ease of use

Nozzle attachment: Easily snaps on and off

Cleaning performance

For best results: Change nozzle every 6 months

Design and finishing

Nozzle color: Grey

