

PHILIPS

sonicare

Interdental nozzles

AirFloss Ultra

2 nozzles

HX8032/23



A triple burst of better gingival health

Efficacy meets compliance for infrequent flossers

For the most effective clean, Philips Sonicare AirFloss Ultra nozzles should be replaced every six months. AirFloss Ultra removes up to 99.9% of plaque biofilm in the treated area.*

Provides a superior clean


- High performance nozzle

Designed around you

- Guidance tip


Highlights

Guidance tip



Helps ensure correct placement. Slide along tooth's surface near the gumline until the tip fits into the groove between the teeth

High performance nozzle



Features a softer tip for better comfort and a wider spray to cover more interproximal tooth and gum surface.

Specifications

Items included

AirFloss Ultra nozzle: 2

Ease of use

Nozzle attachment: Easily snaps on and off

Cleaning performance

For best results: Change nozzle every 6 months

Design and finishing

Nozzle color: Grey

