

PHILIPS

sonicare

Interdental -
Nozzles

AirFloss

w/ 2 Nozzles

HX8012/09



Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

Helps build a healthy habit


- An easy way to improve cleaning between teeth

Easy to use

- Air and micro-droplet technology
- Guidance tip ensures correct placement
- Nozzle designed for easy access to hard-to-reach areas

Highlights

An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Air and micro-droplet tech



Air and micro-droplet technology

Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

Slim, angled nozzle

Slim, angled nozzle designed for easy access to hard-to-reach areas

Specifications

Ease of use

Nozzle attachment: Easily snaps on and off

Items included

AirFloss Nozzle: 2

Cleaning performance

Nozzles: Best results change every 6mos

