

# PHILIPS

## sonicare

Interdental nozzles

AirFloss Ultra

2 nozzles

HX8032/07



## Healthier gums in 2 weeks, guaranteed\*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Ultra nozzles are the easiest way to effectively clean between teeth.

### Improves oral health

- Up to 99.9% plaque removal\*\*\*
- Clinically proven as effective as floss for gum health\*\*
- Helps prevent cavities between teeth

### Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

### The easiest way to effectively clean between teeth

- An easy way to start a healthy habit

# Highlights

## Up to 99.9% plaque removal



AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*\*

## Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

## Better gum health




Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.\*\* Helps improve gum health in as little as two weeks.

## Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

## High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

## Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

# Specifications

## Items included

AirFloss Ultra nozzle: 2

## Ease of use

Nozzle attachment: Easily snaps on and off

## Cleaning performance

For best results: Change nozzle every 6 months

## Design and finishing

Nozzle color: Grey

