

PHILIPS

sonicare

Sonic electric
toothbrush - Trial

FlexCare Platinum

3 modes, 3 intensities

2 brush heads

With pressure sensor



HX9142/31



Their techniques may vary, but their results won't

Our first brush with Adaptive Cleaning Technology

At the core of this technology are two components: the AdaptiveClean brush head and an intuitive pressure sensor. Both features help your patients to brush smarter.

Offers a customized brushing experience

- New Deep Clean: Effortlessly tackle trouble spots

Designed around you

- InterCare brush head

Guides you to follow dental recommendations

- Intuitive pressure sensor

Proven to improve oral health

- Quadpacer
- Smarttimer
- UV sanitizer

Provides a superior clean

- AdaptiveClean brush head
- Customizable brushing experiences
- Patented Sonicare toothbrush technology

Highlights

Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

Deep Clean mode



Our all-new Deep Clean mode delivers an exceptional clean in hard-to-reach places. The unique motion of the brush head in this mode plus an extended three-minute brushing time makes it the ideal choice when you want to pay extra attention to trouble spots.

AdaptiveClean brush head



Soft, flexible sides made of rubber allow the bristles to conform to the unique shape of patients' teeth and gums to provide up to 4x more surface contact*. The brush head's flexibility helps the bristles gently track along the gumline and absorb excessive brushing pressure. Enhanced sweeping motion delivers a unique in-mouth feel and superior cleaning

along the gumline and in hard-to-reach areas*. Proven to remove 10x more plaque than a manual toothbrush in hard-to-reach areas*.

InterCare brush head




Extra-long, high-density bristles reach deeper between teeth to help remove significantly more plaque* interproximally. Ideal for patients with bleeding gums or increased risk of gum inflammation.

Intuitive pressure sensor



If too much pressure is applied, the handle gently pulsates to tell patients they are brushing too hard to guide proper technique. In a study of U.S. patients, 7 out of 10 aggressive brushers tested reduced their brushing pressure when using the pressure sensor.

Custom brushing experiences



Three intensity settings and three cleaning modes provide optimal comfort and performance for a more personalized brushing


experience, encouraging better compliance. The three-minute Deep Clean mode delivers an exceptional clean in hard-to-reach places. It's the perfect mode for patients who need to pay extra attention to trouble spots.

UV sanitizer



Kills up to 99% of bacteria and viruses*. Can be available as an accessory in one of the dispensing SKUs.

Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

Smarttimer



Two-minute timer helps ensure dental professional recommended brushing time.

Specifications

Modes

Clean: For exceptional everyday clean

White: Removes surface stains

Deep Clean: For an invigorating deep clean

3 intensities: High, Low, Medium

Items included

Handles: 1 FlexCare Platinum

Brush heads: 1 AdaptiveClean, 1 InterCare standard

Charger: 1

Travel case: 1

Design and finishing

Color: Platinum gray

Cleaning performance

Speed: Up to 62000 brush movement/min

Performance: Removes up to 10x more plaque*

Health benefits: Improves gum health in only two weeks

Whitening benefits: Helps naturally remove stains

Timer: Quadpacer and SmartTimer

Pressure sensor: Vibrates handle to alert user

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks***

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks**

Battery type: Lithium ION

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty



* Than a DiamondClean brush head

* Than a manual toothbrush

* Data on file

* After a single use of each toothbrush vs. Oral-B 3000/6000 CrossAction

* ** based on two periods of two-minute brushings per day