

# PHILIPS

## sonicare

Sonic electric  
toothbrush with app

DiamondClean  
Smart

Personalized coaching

Smart brush head sensors

Smart brush head recognition

5 modes, 3 intensities



HX9924/06



## Our best ever toothbrush, for complete oral care Feel confident, with 100% coverage

The DiamondClean Smart is our best toothbrush, for complete care. Four high-performance brush heads let you focus on all areas of your oral health, and our Smart Sensor technology gives you personalized feedback and coaching.

### Proven to improve oral health

- Removes up to 10x more plaque\* for a deep clean
- Up to 7x healthier gums\* in just 2 weeks
- Remove up to 100% more stains\*\* in just 3 days
- Enjoy long-lasting fresh breath

### Personalized in-app coaching to perfect brushing

- Guides you to brush and highlights missed areas
- Takes you back to any missed spots for 100% coverage
- Start taking better care of your trouble areas
- Alerts you if you're pressing too hard
- Guides you to reduce your scrubbing

### Optimize your brushing

- Brush heads automatically select the optimal settings
- Make sure you're getting the most from your brush head
- Choose from 5 modes and 3 intensity settings



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## Say goodbye to plaque




Click on the C3 Premium Plaque Control brush head to experience our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact\*\*\* and remove up to 10x more plaque from hard-to-reach areas\*.

## Care for your gums



Click on the G3 Premium Gum Care brush head to improve your gum health. Its smaller size and targeted gum line bristles deliver gentle yet effective cleaning along the gum line, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation\* and up to 7x healthier gums in just two weeks.\*

## Whiter teeth, fast




Click on the W3 Premium White brush head to remove surface stains and reveal your whitest smile. With its densely packed central stain removal bristles, it's clinically proven to remove up to 100% more stains in just three days.\*\*

## TongueCare+ tongue brush




Click on the TongueCare+ tongue brush to gently remove odor causing bacteria from the pores of your tongue. Its 240 specially designed MicroBristles get between all your tongue's ridges and grooves to remove bacteria and debris that lead to bad breath. Team with our antibacterial BreathRx tongue spray for superior cleaning and super fresh breath.

## Never miss a spot




The Philips Sonicare app lets you know when you've achieved a thorough clean, and coaches you to be a better, more mindful brusher.

## TouchUp feature




If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

## Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.

## Too much pressure?



You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your handle will flash. This is a gentle reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.



[asimpleswitch.com](https://www.asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Modes

**Clean:** For exceptional everyday clean

**White+:** To remove surface stains

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

**TongueCare+:** For long lasting clean breath

**3 intensities:** High, Medium, Low

## Items included

**Handle:** 1 DiamondClean Smart

**Brush heads:** 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 W3 Premium White, 1 TongueCare+ tongue brush

**Glass charger:** 1

**Travel case:** USB charging travel case

## BrushSync Mode Pairing

**C3 Premium Plaque Control:** Pairs with Clean mode

**G3 Premium Gum Care brush head:** Pairs with Gum Health mode

**W3 Premium White brush head:** Pairs with White+ mode

**TongueCare+ tongue brush:** Pairs with TongueCare mode

## Connectivity

**Bluetooth® wireless technology:** Connected brushing app

## Compatibility

**Android compatibility:** Android phones, Bluetooth 4.0 enabled tablets

**iOS compatibility:** iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher

## Design and finishing

**Color:** White

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Performance:** Removes up to 10x more plaque\*

**Health benefits:** Up to 7x healthier gums\*\*\*\*

**Whitening benefits:** Up to 100% less stains\*\*\*\*\*

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**TouchUp:** Ensures 100% coverage

**Timer:** BrushPacer and SmarTimer

## Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**Location sensor:** Tracks & improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

**3D mouth map:** Highlights trouble areas

**TouchUp:** To address missed spots

**Position sensor:** Guides to areas you missed

## Ease of use

**Handle compatibility:** Easy click-on brush heads

**Brush head recognition:** Syncs with optimal mode

**Replacement reminder:** To always ensure best results, reminder icon lights up

**Battery indicator:** Illuminated icon indicates battery life

**Handle:** Slim ergonomic design

**Display:** Illuminated display

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** 14 days\*\*\*\*\*

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty



\* Removes up to 7x more plaque vs. a manual toothbrush

\* \* as compared to a manual toothbrush with a leading whitening toothpaste

\* \* compared to DiamondClean

\* \* \* \* in Gum Health mode vs a manual toothbrush at 2 weeks

\* \* \* \* after 3 days of use in White+ mode vs a manual toothbrush

\* \* \* \* based on two periods of two-minute brushings per day on standard mode