



Sonic electric toothbrush with app

- Personalized coaching
- Smart brush head sensors
- Smart brush head recognition
- 5 modes, 3 intensities

HX9924/21

Our best ever toothbrush, for complete oral care

Feel confident, with 100% coverage

The DiamondClean Smart is our best toothbrush, for complete care. Four high-performance brush heads let you focus on all areas of your oral health, and our Smart Sensor technology gives you personalized feedback and coaching.

Benefits

Proven to improve oral health

- Removes up to 10x more plaque* for a deep clean
- Up to 7x healthier gums* in just 2 weeks
- Remove up to 100% more stains** in just 3 days
- Enjoy long-lasting fresh breath

Personalized in-app coaching to perfect brushing

- Start taking better care of your trouble areas
- Guides you to brush and highlights missed areas
- Guides you to reduce your scrubbing

- Alerts you if you're pressing too hard
- Takes you back to any missed spots for 100% coverage


Optimize your brushing

- Brush heads automatically select the optimal settings
- Make sure you're getting the most from your brush head
- Choose from 5 modes and 3 intensity settings
- Easy and stylish charging - home and away

PHILIPS


Features

Say goodbye to plaque



Click on the C3 Premium Plaque Control brush head to experience our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact*** and remove up to 10x more plaque from hard-to-reach areas*.

Care for your gums



Click on the G3 Premium Gum Care brush head to improve your gum health. Its smaller size and targeted gum line bristles deliver gentle yet effective cleaning along the gum line, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation* and up to 7x healthier gums in just two weeks.*

Whiter teeth, fast



Click on the W3 Premium White brush head to remove surface stains and reveal your whitest smile. With its densely packed central stain removal bristles, it's clinically proven to remove up to 100% more stains in just three days.**

TongueCare+ tongue brush




Click on the TongueCare+ tongue brush to gently remove odor causing bacteria from the pores of your tongue. Its 240 specially designed MicroBristles get between all your tongue's ridges and grooves to remove bacteria and debris that lead to bad breath. Team with our antibacterial BreathRx tongue spray for superior cleaning and super fresh breath.

Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.

Never miss a spot



The Philips Sonicare app lets you know when you've achieved a thorough clean, and coaches you to be a better, more mindful brusher.

Easy does it



With an electric toothbrush, you let the brush do the work. We've built a scrubbing sensor into your handle as a gentle reminder to reduce scrubbing. This way, you can improve your technique and get a gentler, more effective clean.

BrushSync mode pairing



Wondering which mode and intensity you should use? Wonder no longer. A microchip in your brush head tells your DiamondClean Smart which one you're using. So if you click on a gum-care brush head, for example, your toothbrush knows to select the optimal mode and intensity to provide a gentle yet effective care for your gums. All you need to do is to press the power button.

Smart replacement reminders



Brush heads wear-out over time and become less effective at removing plaque. With smart brush head recognition you'll never lose track of when to replace your brush head. Your toothbrush tracks how long and hard you've brushed with each brush head, and will alert you when it's time to replace. You can also monitor the lifetime of your brush heads in the Philips Sonicare app, and even order replacements so you're never caught short.

5 modes, 3 intensity settings




DiamondClean Smart comes with five modes which cover a range of brushing needs: Clean mode is for exceptional daily cleaning, White+ is for removing stains, Deep Clean+ gives you an invigorating deep clean, Gum Health provides a gentle yet effective clean for gums, and Tongue Care is the perfect

DiamondClean Smart 9500


mode for cleaning your tongue. It breaks up tongue debris and sweeps away bad-breath bacteria. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

Too much pressure?



You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your handle will flash. This is a gentle reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

TouchUp feature



If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

Design with storage in mind



You'll enjoy two weeks of regular use from a single full charge. However, if you run low while traveling, your deluxe USB travel case doubles as a charger. Place your toothbrush in the case and either connect it to your laptop or a wall socket. As for charging at home - simply place your brush in your sleek charger glass - it fits stylishly into your bathroom and also doubles as a rinsing glass.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Modes

Clean	For exceptional everyday clean
White+	To remove surface stains
Deep Clean+	For an invigorating deep clean
Gum Health	Special attention to molars
TongueCare+	For long lasting clean breath
3 intensities	High Medium Low

iOS compatibility

Bluetooth 4.0 enabled tablets
iPad 3rd Gen or higher
iPhone 4S or higher
with iOS7 or higher

Design and finishing

Color	Pink
-------	------

Cleaning performance

Speed	Up to 62000 brush movement/min
Performance	Removes up to 10x more plaque*
Health benefits	Up to 7x healthier gums****
Whitening benefits	Up to 100% less stains*****
Pressure feedback	Ring lights up in purple Vibrates handle to alert user
TouchUp	Ensures 100% coverage
Timer	BrushPacer and SmarTimer

Smart sensor technology

Pressure sensor	Alerts when brushing too hard
Location sensor	Tracks & improves coverage
Scrubbing sensor	Guides to reduce scrubbing
3D mouth map	Highlights trouble areas
TouchUp	To address missed spots
Position sensor	Guides to areas you missed

Ease of use

Handle compatibility	Easy click-on brush heads
Brush head recognition	Syncs with optimal mode
Replacement reminder	To always ensure best results reminder icon lights up
Battery indicator	Illuminated icon indicates battery life
Handle	Slim ergonomic design
Display	Illuminated display

Items included

Handle	1 DiamondClean Smart
Brush heads	1 C3 Premium Plaque Control 1 G3 Premium Gum Care 1 W3 Premium White 1 TongueCare+ tongue brush
Glass charger	1
Travel case	USB charging travel case

BrushSync Mode Pairing

C3 Premium Plaque Control	Pairs with Clean mode
G3 Premium Gum Care brush head	Pairs with Gum Health mode
W3 Premium White brush head	Pairs with White+ mode
TongueCare+ tongue brush	Pairs with TongueCare mode

Connectivity

Bluetooth® wireless technology	Connected brushing app
--------------------------------	------------------------

Compatibility

Android compatibility	Android phones
-----------------------	----------------

DiamondClean Smart 9500

Technical specifications

Battery	Rechargeable
Battery type	Lithium ION
Operating time (full to empty)	14 days*****

Power

Voltage	110-220 V
---------	-----------

Service

Warranty	2-year limited warranty
----------	-------------------------

* Removes up to 7x more plaque vs. a manual toothbrush
* * as compared to a manual toothbrush with a leading whitening toothpaste
* * compared to DiamondClean
* *** in Gum Health mode vs a manual toothbrush at 2 weeks
* **** after 3 days of use in White+ mode vs a manual toothbrush
* **** based on two periods of two-minute brushings per day on standard mode



data subject to change
2021, September 30
Version: 15.5.1
UPC: 00075020077080

© 2021 Koninklijke Philips N.V.
All rights reserved.
Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.
www.philips.com