

# PHILIPS

## sonicare



### Oral Irrigator System

#### Power Flosser 7000 System

Quad Stream technology

Whisper quiet design

Pulse waves guide as you go  
with ExpertClean 7300



HX3921/40

## Up to 9X as effective in removing plaque\*

### An effortless, thorough clean

Enjoy an effortless, thorough clean between your teeth and improve your gum health. Our Quad Stream technology helps you cover more area with less effort, while Pulse Wave technology gently guides you from tooth to tooth. Easy.

#### **An effortless clean with advanced technology**

- Easy to fill 600ml reservoir
- Quadstream technology: removes plaque from below the gumline
- Quad Stream technology for faster, more effective flossing
- Standard nozzle for a regular clean
- Comfort nozzle for sensitive gums
- Pulse Wave technology guides from tooth to tooth
- Whisper quiet for great results without disturbing others
- Quad Stream nozzle for maximum coverage
- Easy magnetic dock with auto shut-off
- Compact nozzle case
- Philips Sonicare ExpertClean 7300

#### **Customize easily to your needs**

- 4 flossing modes, 10 intensities


# Power Flosser 7000 System

Oral Irrigator System

HX3921/40


## Highlights

### Covers more area, effortlessly




X-shaped water streams cover more surface area for a faster, deeper clean versus manual floss. Effortlessly.

### Cleans from below the gumline



The Quadstream technology removes up to 99% of plaque from 6mm deep periodontal pockets.

### Guides from tooth to tooth



Gentle pulses of water guide you from tooth-to-tooth so that you don't miss a spot.

### Whisper quiet design



Engineered for a discreet clean, so you get great results without disturbing others.

### 4 modes, 10 intensities



Enjoy a clean that's tailored to your needs. Clean mode uses a continuous flow of water for a great everyday clean; Deep Clean uses distinct pulses for a more thorough clean; Sensitive mode is extra gentle for sensitive gums; Massage mode soothes and stimulates the gumline. Whenever you need, adjust the intensity for your comfort using the LED light display as your guide.


### F3 Quad Stream nozzle



For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams


between teeth and along the gum line. The soft rubber guide is gentle on gums and every nozzle simply clicks on and off.

### F2 Comfort nozzle



Want to take care of sensitive gums? Click on the Comfort nozzle and gently guide the soft rubber tip along your gumline.

### F1 Standard nozzle



Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

### Magnetic dock w/ auto shut-off



When you're finished, simply return the handle to the easy magnetic dock. The handle automatically switches from pause to off.

# Power Flosser 7000 System

Oral Irrigator System

HX3921/40

## Highlights

### 600ml resevoir



The 600ml resevoir with easy flip lid holds enough water for one cleaning. To refill simply lift out, fill, and click back in.

### Compact nozzle case

Store extra nozzles hygienically in our compact and stylish nozzle case.

### ExpertClean 7300



Complete your clean with the Philips Sonicare ExpertClean 7300. It's designed to give you an exceptional daily clean and help you keep up

good brushing habits. Smart sensors help you maintain the perfect pressure and an in-app progress report encourages you to keep up your routine. When you're done, simply return the handle to the built-in charger, or use the travel charger when you're on the go.

# Power Flosser 7000 System

Oral Irrigator System

HX3921/40

## Specifications

### Power

**Voltage:** Multi-voltage charger

### Service

**Warranty:** 2-year limited warranty

### Ease of use

**Handle:** Slim ergonomic design

**Nozzle attachment:** Easily clicks on and off

### Performance

**Cleaning:** Clean entire mouth in 60-90sec

### Items included

**Handle:** 1 ExpertClean

**Travel case:** 1

**Charger:** 1

**Brush heads:** 1 C3 Premium Plaque Control

**F1 Standard nozzle:** 1

**F3 Quad Stream nozzle:** 1

**F2 Comfort nozzle:** 1

**Power Flosser:** 1

**Nozzle storage case:** 1

### Modes

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Massage:** Invigorating massage

**Intensities:** 10

**Sensitive:** Gentle teeth and gum cleaning

### Health benefits

**Gum health:** Up to 270% more effective\*\*

**Plaque removal:** 9x as effective\*

**Sub-gingival:** Removes up to 99% of plaque from below the gumline\*\*\*



\* in hard-to-reach areas vs. manual toothbrush + floss. Using Sonicare Power Flosser with Quad Stream tip on setting 8 and Sonicare ExpertClean power toothbrush with Premium Plaque Control brush head in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss

\* \*vs. manual toothbrush + floss. Using Sonicare Power Flosser with Quad Stream tip on setting 8 and Sonicare ExpertClean power toothbrush with Premium Plaque Control brush head in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss

\* \*\*In a lab study (actual results may vary), in the deepest 2mm of a 6mm periodontal pocket model.