

# PHILIPS

## Power Toothbrush

### Philips One by Sonicare

The One that cleans better

The One for on-the-go

The One without hassle

The One in multiple colors



HY1200/06



## One up your brushing

The new electric One cleans your teeth better for a brighter, whiter smile. Wherever you go, this One goes with you. Its sleek, lightweight design and portable case makes it the perfect brush for you. Which color is your One?

#### Better brushing

- Good vibrations

#### On-the-go

- Travel case to take with you wherever whenever

#### Hassle-free

- Routine made easy with brushing guidance
- One charge, one month of brushing
- Keep your brush fresh

#### Multiple colors

- A color for every style

# Highlights

## Good vibrations



Think of the One like a manual toothbrush that works harder and faster. 13,000 microvibrations per minute and contoured bristles gently clean your teeth for a whiter, brighter smile.

## The One that goes with you



The One is made to travel so you can get that clean feeling any place. It's sleek and lightweight and fits neatly into a compact travel case.

## Routine made easy




Sometimes, we all need a little guidance. The One buzzes every 30 seconds encouraging you to move on to the next part of your mouth. A longer buzz at 2 minutes lets you know when the job is done.

## A color for every style



From Shadow to Shimmer, Snow to Sage, choose the color to match your style. You can mix and match the brush head with the handle and create your own color combination.

## One charge, one month of brush



Charging your One is easy. Just pop off the end cap, plug in your USB-C cable and wait until the light indicator stops flashing. With a single charge, you can brush for 30 days. That's one month of whiter, brighter smiling.

## Keep it fresh



Dentists recommend changing your brush head every three months for the best results. Philips One brush heads come in a variety of vibrant colors, so you can update your look as well as your brush head every time.

# Specifications

## Technical specifications

Operating time (full to empty): Up to 30 days

## Design and finishing

Color: Shadow

## Compatibility

**Brush head compatibility:** Only compatible with, Philips One handles

**Handle compatibility:** Not compatible with other, Philips Sonicare brush heads

## Items included

**Travel case:** Philips One travel case

**Handle:** 1 Philips One

**Brush heads:** 1 Philips One brush head

