

PHILIPS

sonicare

Tongue care kit

BreathRx

2 oz tongue spray

2 tongue cleaners

DIS359/03



Attack and neutralize bad breath at the source Part of the #1 recommended breath care system

A cleaner tongue for fresher breath. The Philips Sonicare BreathRx tongue spray destroys bad breath bacteria and neutralizes odors, while the tongue cleaner gently removes odor-causing bacteria build-up from the tongue.

A custom clean

- Soft and regular side

A powerful combination of ingredients

- Kills 99% of bad breath bacteria
- Neutralizes odors

The 3 step system

- Lasting fresh breath in 3 easy steps
- Step 2: Scrape

Highlights

Step 2




Up to 80% of bad breath comes from odor-causing bacteria that hides in the uneven surface of the tongue. The active ingredient in the BreathRx tongue spray kills bad breath bacteria on contact and neutralizes odor-causing gases. The flexible BreathRx tongue cleaner gently removes bacteria build-up from the tongue.

2 Sides



Two cleaning sides for a customizable cleaning experience. Choose Regular for a deep clean or Soft for sensitive tongues.

CPC



CPC, a powerful antibacterial ingredient, attacks and kills bad breath bacteria on contact.

Zytex



Zytex, our unique blend of zinc, and thymol & eucalyptus essential oils, neutralizes bad breath odors for instantly fresher breath.

3 Steps



Don't just mask odors, eliminate them in 3 easy steps - brush, scrape, rinse. Step 1: Brush away odor-causing plaque with BreathRx Whitening toothpaste. Step 2: Get to the bacteria that hides within the uneven surface of your tongue with the BreathRx tongue spray and tongue cleaner. The active ingredient in the spray kills bad breath bacteria on contact and neutralizes odors. The flexible tongue cleaner gently removes bacteria build-up from the tongue. Step 3: Rinse away the remaining germs and debris with the BreathRx antibacterial mouth rinse.

Specifications

Items included

Antibacterial tongue spray: 2oz
Tongue cleaners: 2

Ingredients

Antibacterial tongue spray: 0.09% CPC, Zytex

Fresh breath system

Step: 2

Expiration

After opened: 6 months

Shelf life: 36 months

