



Philips Sonicare AirFloss
Interdental - Nozzles

w/ 2 Nozzles

HX8012/05

Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

Easy to use

- Air and micro-droplet technology
- Guidance tip ensures correct placement
- Nozzle designed for easy access to hard-to-reach areas

Helps build a healthy habit

- An easy way to improve cleaning between teeth

PHILIPS
sonicare

Specifications

Ease of use

- Nozzle attachment: Easily snaps on and off

Items included


- AirFloss Nozzle: 2

Cleaning performance

- Nozzles: Best results change every 6mos


Highlights

An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Air and micro-droplet tech



AirFloss utilizes unique air and micro-droplet technology to send a burst of pressurized air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

Slim, angled nozzle

Slim, angled nozzle designed for easy access to hard-to-reach areas



Issue date 2022-08-16

Version: 5.5.2

EAN: 08 71010 36240 11

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* than manual toothbrushing alone