

# PHILIPS

## sonicare

Sonic electric  
toothbrush with app

DiamondClean  
Smart

Personalized coaching

Smart brush head sensors

Smart brush head recognition

4 modes, 3 intensities



HX9902/65



## Our best ever toothbrush, for complete oral care Feel confident, with 100% coverage

The DiamondClean Smart is our best toothbrush, for complete care. High-performance brush heads let you focus on all areas of your oral health, and our Smart Sensor technology gives you personalized feedback and coaching.

### Proven to improve oral health

- Removes up to 10x more plaque\* for a deep clean

### Personalized in-app coaching to perfect brushing

- Guides you to brush and highlights missed areas
- Takes you back to any missed spots for 100% coverage
- Start taking better care of your trouble areas
- Alerts you if you're pressing too hard
- Guides you to reduce your scrubbing

### Optimize your brushing

- Easy and stylish charging - home and away
- Brush heads automatically select the optimal settings
- Choose from 4 modes and 3 intensity settings
- Make sure you're getting the most from your brush head

### ExpertClean users experience improved oral health

- Up to 7x healthier gums\* in just 2 weeks


# Highlights

## Say goodbye to plaque




Click on the C3 Premium Plaque Control brush head to experience our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact\*\*\* and remove up to 10x more plaque from hard-to-reach areas\*.

## Up to 7x healthier gums\*




Focus on improving your gum health with the G3 Premium Gum Care brush head. Its smaller size and targeted gum line bristles make it the perfect way to keep this area gum disease free. People who use this brush head have up to 100% less gum inflammation\* and up to 7x healthier gums in just two weeks.\*

## 4 modes, 3 intensity settings




With DiamondClean Smart you get the full works. Our 3 intensity settings let you boost your clean, and our 4 modes cover all your brushing needs; Clean mode for exceptional daily cleaning, White+ to remove stains, Deep Clean+ for an invigorating deep clean, and Gum Health to gently clean along your gum line.

## Never miss a spot




The Philips Sonicare app lets you know when you've achieved a thorough clean, and coaches you to be a better, more mindful brusher.

## TouchUp feature




If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

## Smart replacement reminders




Brush heads wear-out over time and become less effective at removing plaque. With smart brush head recognition you'll never lose track of when to replace your brush head. Your toothbrush tracks how long and hard you've brushed with each brush head, and will alert you when it's time to replace. You can also monitor the lifetime of your brush heads in the Philips Sonicare app, and even order replacements so you're never caught short.

## Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.

## Too much pressure?



You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your handle will flash. This is a gentle reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

## Easy does it



With an electric toothbrush, you let the brush do the work. We've built a scrubbing sensor into your handle as a gentle reminder to reduce scrubbing. This way, you can improve your technique and get a gentler, more effective clean.

# Specifications

**Connectivity**

**Bluetooth® wireless technology:** Connected brushing app

**Power**

**Voltage:** 110–220 V

**Technical specifications**

**Operating time (full to empty):** 14 days\*\*\*\*\*

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Design and finishing**

**Color:** Pastel pink

**Service**

**Warranty:** 2-year limited warranty

**Compatibility**

**Android compatibility:** Android phones, Bluetooth 4.0 enabled tablets

**iOS compatibility:** iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher

**Ease of use**

**Handle:** Slim ergonomic design

**Battery indicator:** Illuminated icon indicates battery life

**Handle compatibility:** Easy click-on brush heads

**Display:** Illuminated display

**Brush head recognition:** Syncs with optimal mode

**Replacement reminder:** To always ensure best results, reminder icon lights up

**Items included**

**Brush heads:** 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care

**Travel case:** USB charging travel case

**Handle:** 1 DiamondClean Smart

**Glass charger:** 1

**Cleaning performance**

**Health benefits:** Helps improve gum health

**Performance:** Removes up to 10x more plaque\*

**Speed:** Up to 62000 brush movement/min

**Timer:** BrushPacer and SmarTimer

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**TouchUp:** Ensures 100% coverage

**Modes**

**3 intensities:** Low, Medium, High

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

**White+:** To remove surface stains

**BrushSync Mode Pairing**

**C3 Premium Plaque Control:** Pairs with Clean mode

**G3 Premium Gum Care brush head:** Pairs with Gum Health mode

**W3 Premium White brush head:** Pairs with White+ mode

**Smart sensor technology**

**Pressure sensor:** Alerts when brushing too hard

**Location sensor:** Tracks & improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

**3D mouth map:** Highlights trouble areas

**TouchUp:** To address missed spots



\* Removes up to 7x more plaque vs. a manual toothbrush  
\* \* as compared to a manual toothbrush with a leading whitening toothpaste  
\* \* compared to DiamondClean  
\* \*\*\* in Gum Health mode vs a manual toothbrush at 2 weeks  
\* \*\*\*\* after 3 days of use in White+ mode vs a manual toothbrush  
\* \*\*\*\* based on two periods of two-minute brushings per day on standard mode