

# PHILIPS

## sonicare



Oral Irrigator

Power Flosser 5000

Quad Stream technology

Whisper quiet design

Pulse waves guide as you go



HX3811/20

## Up to 180% more effective for healthier gums\*

### An effortless, thorough clean

Enjoy an effortless, thorough clean between your teeth and improve your gum health. Our Quad Stream technology helps you cover more area with less effort, while Pulse Wave technology gently guides you from tooth to tooth. Easy.

#### **An effortless clean with advanced technology**

- Easy to fill 550ml reservoir
- Quadstream technology: removes plaque from below the gumline
- Quad Stream technology for faster, more effective flossing
- Standard nozzle for a regular clean
- Pulse Wave technology guides from tooth to tooth
- Quad Stream nozzle for maximum coverage
- Easy magnetic dock with auto shut-off
- Whisper quiet for great results without disturbing others

#### **Customize easily to your needs**

- 2 flossing modes, 10 intensities


# Power Flosser 5000

Oral Irrigator

HX3811/20

## Highlights

### Covers more area, effortlessly




X-shaped water streams cover more surface area for a faster, deeper clean versus manual floss. Effortlessly.

### Cleans from below the gumline



The Quadstream technology removes up to 99% of plaque from 6mm deep periodontal pockets.

### Guides from tooth to tooth



Gentle pulses of water guide you from tooth-to-tooth so that you don't miss a spot.

### Whisper quiet design


Engineered for a discreet clean, so you get great results without disturbing others.

### 2 modes, 10 intensities



Enjoy a clean that's tailored to your needs. Clean mode uses a continuous flow of water for a great everyday clean, while Deep Clean uses distinct pulses for a more thorough clean. Whenever you need, adjust the intensity to your ideal comfort level.

### F1 Standard nozzle



Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

### F3 Quad Stream nozzle



For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams

between teeth and along the gum line. The soft rubber guide is gentle on gums and every nozzle simply clicks on and off.

### 550ml reservoir



The 550ml reservoir holds enough water for one cleaning. To refill simply lift out, fill and click back in.

### Magnetic dock w/ auto shut-off



When you're finished, simply return the handle to the easy magnetic dock. The handle automatically switches from pause to off.

# Power Flosser 5000

Oral Irrigator

HX3811/20

## Specifications

### Power

**Voltage:** Multi-voltage charger

### Service

**Warranty:** 2-year limited warranty

### Ease of use

**Handle:** Slim ergonomic design

**Nozzle attachment:** Easily clicks on and off

### Performance

**Cleaning:** Clean entire mouth in 60-90sec

### Items included

**F1 Standard nozzle:** 1

**F3 Quad Stream nozzle:** 1

**Power Flosser:** 1

### Modes

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Intensities:** 10

### Health benefits

**Gum health:** Up to 180% more effective\*

**Plaque removal:** Up to 3X as effective\*\*,

Removes up to 99.9% plaque\*\*\*

**Sub-gingival:** Removes up to 99% of plaque from below the gumline\*\*\*\*

© 2025 Koninklijke Philips N.V.  
All Rights reserved.

Issue date 2025-04-21  
Version: 5.5.3

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

UPC: 0 75020 09115 4

[www.philips.com](http://www.philips.com)



\* vs. floss. When used with Quad Stream nozzle on setting 8 with a manual toothbrush in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss

\* \*when used with Quad Stream tip on setting 8 with a manual toothbrush in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss.

\* \* \*in a lab study. Actual in-mouth results may vary.

\* \* \* \*In a lab study (actual results may vary), in the deepest 2mm of a 6mm periodontal pocket model.