

# PHILIPS



Sonic electric  
toothbrush with app

ExpertClean 7300

Connected brushing made easy

Built-in pressure sensor

Smart brush head recognition

3 modes, 3 intensities



HX9610/17

## Everything you need for great oral health

Up to 7x healthier gums\*

Keeping up good brushing habits is tough. ExpertClean will guide you between dentist checkups. Know when to adjust your pressure with built-in smart sensors and use the App's Progress Report to help you stick to your brushing routine.

### Enjoy a stylish and user-friendly design

- Simple to charge and style to boot
- Makes traveling easier for you

### ExpertClean users experience improved oral health

- Hard on plaque, gentle on your gums
- Deepest clean while protecting your sensitive gums
- Removes up to 10x more plaque\* for a deep clean
- Up to 7x healthier gums\* in just 2 weeks

### Let ExpertClean adapt to your needs

- Personalized brushing experience
- Lets you know if you're pressing too hard
- Always know when to replace your brush heads
- BrushSync automatically selects the best mode for you

### Stay on track with personalized guidance

- Guidance to start and keep up healthy habits

# ExpertClean 7300

Sonic electric toothbrush with app

HX9610/17


## Highlights

### Removes up to 10x more plaque\*




Know you're getting the deepest possible clean with The C3 Premium Plaque Control brush head. The soft, flexible bristles are designed to curve around the contours of your teeth, giving you 4x more surface contact\*\* and up to 10x more plaque removal from hard-to-reach spots.

### Up to 7x healthier gums\*



Focus on improving your gum health with the G3 Premium Gum Care brush head. Its smaller size and targeted gum line bristles make it the perfect way to keep this area gum disease free. People who use this brush head have up to 100% less gum inflammation\* and up to 7x healthier gums in just two weeks.\*

### Progress Report



ExpertClean gives you the guidance you need to improve and maintain healthy brushing habits in between dentist checkups. Built-in smart sensors let you know when you're using too much pressure, and by connecting your brushing

experience to the Sonicare app, a personalized Progress Report helps you stay on track, to see how much you've improved over time.

### 3 modes, 3 intensities




ExpertClean comes with Clean, Gum Health and Deep Clean+ modes to take care of your brushing needs. Clean mode is for exceptional everyday cleaning, Gum Health provides a gentle yet effective clean for gums and Deep Clean+ gives you an invigorating deep clean. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

### Smart brush head recognition



Smart brush heads ensure you're using the right mode and intensity for the best possible clean. For an example, say you're using the G3 Premium Gum Care brush head. Your ExpertClean's BrushSync technology will automatically sync your brush head with the Gum Health mode to help give your gums a clean bill of health.

### Built-in pressure sensor



You might not notice if you're brushing too hard, but your ExpertClean will. If you need to ease up, this intelligent toothbrush will make a pulsing sound. It's a heads up to let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

### BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

### C3 Premium Plaque Control



The C3 Premium Plaque Control brush head is designed to give you your deepest clean yet. Soft flexible sides and bristles perfectly curve around the contours of your teeth, giving you 4x more surface contact and helping target hard-to-reach spots\*.

# ExpertClean 7300

Sonic electric toothbrush with app

HX9610/17

## Highlights

### G3 Premium Gum Care brush head



The G3 Premium Gum Care brush head is designed to help you improve your gum health. Its smaller size and targeted gumline bristles gently and effectively clean the gum line, where gum disease commonly starts out.

### Long-life battery



Your ExpertClean will last for up to two weeks on a single charge with regular use.

### Makes traveling easier for you



Our premium travel case lets you store your toothbrush hygienically, while our compact charging base keeps you topped up when you're on the go. You do get two weeks of regular use from a single full charge, but the charger's just the ticket for longer trips.

# ExpertClean 7300

Sonic electric toothbrush with app

HX9610/17

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Voltage:** 110-220 V

### Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** 14 days\*\*\*

**Battery type:** Lithium ION

### Design and finishing

**Color:** Black

### Service

**Warranty:** 2-year limited warranty

### Compatibility

**Android compatibility:** Android phones, Bluetooth 4.0 enabled tablets

**iOS compatibility:** iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

### Ease of use

**Handle:** Slim ergonomic design

**Handle compatibility:** Easy click-on brush heads

**Replacement reminder:** To always ensure best results, reminder icon lights up

### Items included

**Handle:** 1 ExpertClean

**Travel case:** 1

**Charger:** 1

**Brush heads:** 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care

### Cleaning performance

**Performance:** Removes up to 10x more plaque\*

**Health benefits:** Up to 7x healthier gums\*

**Timer:** BrushPacer and SmarTimer

**Pressure feedback:** Vibrates handle to alert user

### Modes

**3 intensities:** Low, Medium, High

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

### BrushSync Mode Pairing

**C3 Premium Plaque Control:** Pairs with Clean mode

**G3 Premium Gum Care brush head:** Pairs with Gum Health mode

### Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**BrushSync Replacement Reminder:** Always know when to, replace brush heads

**BrushSync Technology:** Connects smart handle and, smart brush head

**Location sensor:** Tracks & improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

