

Philips Sonicare Power
Flosser 7000
Oral Irrigator

Quad Stream technology
4 flossing modes, 10 intensities
Whisper quiet design
Pulse waves guide as you go



HX3911/30

Flossing, reinvented. With Quad Stream technology

Up to 180% more effective for healthier gums*

Philips Sonicare Power Flosser is flossing, reinvented. The X-shaped stream of water reaches deep between teeth for a gentle, thorough clean. Amazingly clean, amazingly easy.

An effortless clean with advanced technology

- Easy to fill 600ml reservoir
- Quad Stream technology for faster, more effective flossing
- Standard nozzle for a regular clean
- Pulse Wave technology guides from tooth to tooth
- Comfort nozzle for sensitive gums
- Whisper quiet for great results without disturbing others
- Quad Stream nozzle for maximum coverage
- Easy magnetic dock with auto shut-off
- Compact nozzle case

Customize easily to your needs

- 4 flossing modes, 10 intensities

PHILIPS


sonicare

Oral Irrigator

Quad Stream technology 4 flossing modes, 10 intensities, Whisper quiet design, Pulse waves guide as you go

HX3911/30

Highlights

Covers more area, effortlessly



X-shaped water streams cover more surface area for a faster, deeper clean versus manual floss. Effortlessly.

4 modes, 10 intensities




Enjoy a clean that's tailored to your needs. Clean mode uses a continuous flow of water for a great everyday clean; Deep Clean uses distinct pulses for a more thorough clean; Sensitive mode is extra gentle for sensitive gums; Massage mode soothes and stimulates the gumline. Whenever you need, adjust the intensity for your comfort using the LED light display as your guide.

Whisper quiet design



Engineered for a discreet clean, so you get great results without disturbing others.

Guides from tooth to tooth




Gentle pulses of water guide you from tooth-to-tooth so that you don't miss a spot.

F3 Quad Stream nozzle




For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams between teeth and along the gum line. The soft rubber guide is gentle on gums and every nozzle simply clicks on and off.

F2 Comfort nozzle



Want to take care of sensitive gums? Click on the Comfort nozzle and gently guide the soft rubber tip along your gumline.

F1 Standard nozzle



Use the standard nozzle for a regular clean. It uses a single stream of water to remove debris and clean thoroughly between your teeth.

Magnetic dock w/auto shut-off



When you're finished, simply return the handle to the easy magnetic dock. The handle automatically switches from pause to off.

600ml reservoir



The 600ml reservoir with easy flip lid holds enough water for one cleaning. To refill simply lift out, fill, and click back in.

Compact nozzle case

Store extra nozzles hygienically in our compact and stylish nozzle case.

Specifications

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

Performance

- Cleaning: Clean entire mouth in 60-90sec

Items included

- F1 Standard nozzle: 1
- F3 Quad Stream nozzle: 1

- F2 Comfort nozzle: 1
- Power Flosser: 1
- Nozzle storage case: 1

Modes

- Clean: For exceptional everyday clean
- Deep Clean+: For an invigorating deep clean
- Intensities: 10
- Massage: Invigorating massage
- Sensitive: Gentle teeth and gum cleaning

Health benefits

- Gum health: Up to 180% more effective*
- Plaque removal: Up to 3X as effective**, Removes up to 99.9% plaque***



Issue date 2024-08-21

Version: 6.6.3

UPC: 0 75020 09116 1

© 2024 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* vs. floss. When used with Quad Stream nozzle on setting 8 with a manual toothbrush in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss

* *when used with Quad Stream tip on setting 8 with a manual toothbrush in patients with moderate to severe gingivitis vs. using manual toothbrush with string floss.

* **in a lab study. Actual in-mouth results may vary.