

Philips Sonicare AirFloss  
Ultra  
Interdenal nozzles

2 nozzles

HX8032/23

## A triple burst of better gingival health

Efficacy meets compliance for infrequent flossers

For the most effective clean, Philips Sonicare AirFloss Ultra nozzles should be replaced every six months. AirFloss Ultra removes up to 99.9% of plaque biofilm in the treated area.\*

**Designed around you**

- Guidance tip

**Provides a superior clean**

- High performance nozzle

**PHILIPS**

**sonicare**

Interdental nozzles  
2 nozzles

# Specifications

## Design and finishing

- Nozzle color: Grey

## Ease of use

- Nozzle attachment: Easily snaps on and off

## Items included

- AirFloss Ultra nozzle: 2


## Cleaning performance

- For best results: Change nozzle every 6 months

HX8032/23

# Highlights

## Guidance tip



Helps ensure correct placement. Slide along tooth's surface near the gumline until the tip fits into the groove between the teeth

## High performance nozzle



Features a softer tip for better comfort and a wider spray to cover more interproximal tooth and gum surface.



Issue date 2022-07-02

Version: 2.2.1

EAN: 08 71010 38361 00

© 2022 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* In a lab study, actual in-mouth results may vary.