

PHILIPS

sonicare



Bluetooth® connected
toothbrush-Dispense

FlexCare Platinum
Connected

Built-in Bluetooth®

Smart brushing sensors

3 modes, 3 intensities

3 brush heads & UV Sanitizer



HX9193/04

Reveal their habits for improved brushing routines App turns brushing data into feedback and coaching

The new Philips Sonicare FlexCare Platinum Connected, powered by Smart Sensor Technology, reveals your patients' brushing habits and helps coach them into a better oral care routine.

Provides a superior clean

- AdaptiveClean brush head
- BrushPacer and SmarTimer for thorough brushing
- Customizable modes and intensities
- Extended battery life and convenient travel case
- InterCare brush head
- Exclusive bonus items only available to dental professionals
- UV sanitizer

Track your brushing with real-time guidance

- State-of-the-art 3D mouth-map
- Brush head monitor
- Goal-setting and progress reports
- Powered by Smart Sensor Technology
- Highlights trouble spots
- Custom guidance with Touch-Up

FlexCare Platinum Connected

Bluetooth® connected toothbrush-Dispense

HX9193/04

Highlights

Bonus items not sold at retail



Now you can offer these exclusive Philips Sonicare Flexcare Platinum Connected bonus items your patients won't find anywhere else. Designed around patient needs for dental offices like yours, the package includes: 1 ProResults gum health brush head and 1 trial size (2 oz.) BreathRx rinse.

Smart Sensor Technology



Three smart sensors send personalized brushing data to the Philips Sonicare app. Location sensor shows patients where they're brushing too little for better coverage. Scrubbing sensor helps guide patients to reduce scrubbing for optimized technique and a deeper, more gentle clean. Pressure sensor tracks brushing pressure and gently pulsates the handle to tell patients when to adopt a softer touch.

State-of-the-art 3D mouth-map




A personalized 3D mouth-map provides post-brushing analysis so patients can track and refine their routine and identify trouble spots that require extra attention.

Touch-Up



The Touch-Up feature displays customized guidance based on actual brushing data in order to immediately address missed areas.

Highlights trouble spots



Based on advice from their dental professional, patients can highlight and address areas on the personalized mouth-map in need of extra care due to plaque buildup, gum recession or another issue.

Brush head monitor




Using actual brushing time and pressure, the brush head monitor alerts patients in three months or less when it's time to replace their brush head, ensuring continued effectiveness for the best results possible.

Goal-setting, progress reports



The Philips Sonicare app maintains a running history of a patient's personalized brushing data, allowing them to easily review performance and set daily and long-term goals. After 30 days, patients can generate an in-app progress report summarizing their brushing performance, with the added option of sharing reports with their dental professional via email.

BrushPacer and SmarTimer



The BrushPacer tracks brushing through six segments of the mouth for a complete clean, while the SmarTimer ensures that patients brush for the full two-minutes as recommended by dental professionals.


FlexCare Platinum Connected

Bluetooth® connected toothbrush-Dispense

HX9193/04

Highlights

AdaptiveClean brush head



Soft, flexible sides made of rubber allow the bristles to conform to the unique shape of patients' teeth and gums to provide up to 4x more surface contact*. The brush head's flexibility helps the bristles gently track along the gumline and absorb excessive brushing pressure. Enhanced sweeping motion delivers a unique in-mouth feel and superior cleaning along the gumline and in hard-to-reach areas*. Proven to remove 10x more plaque than a manual toothbrush in hard-to-reach areas*.

InterCare brush head



Extra-long, high-density bristles reach deeper between teeth to help remove significantly more plaque* interproximally. Ideal for patients with bleeding gums or increased risk of gum inflammation.

Customizable modes



With three modes, Clean, White and Deep Clean, and three intensity settings, Philips Sonicare FlexCare Platinum Connected provides a customized experience for optimal comfort.

Battery life and travel case



When fully charged, Philips Sonicare FlexCare Platinum Connected retains enough battery life for up to three weeks of use. And with the travel case and charger, on-the-go brushing is easy.

UV sanitizer



Kills up to 99% of bacteria and viruses (E. coli, HSV1, H.A.)

FlexCare Platinum Connected

Bluetooth® connected toothbrush-Dispense

HX9193/04

Specifications

Connectivity

Bluetooth® wireless technology: Connected brushing app

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: White

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 operational system

Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Display: Illuminated display

Brushing time: Up to 2 weeks

Items included

Handles: 1 FlexCarePlatinum Connected

Brush heads: 1 AdaptiveClean, 1 InterCare standard, 1 ProResults gum health

Travel case: 1

Charger: 1

UV sanitizer: 1

Mouth rinse: 2oz BreathRx

Cleaning performance

Health benefits: Improves gum health in only two weeks

Performance: Remove up to 10x more plaque**

Pressure sensor: Vibrates handle to alert user

Whitening benefits: Helps naturally remove stains

Timer: BrushPacer and SmarTimer

Speed: Up to 62000 brush movement/min

Modes

Clean: For exceptional everyday clean

3 intensities: Low, Medium, High

Deep Clean: For an invigorating deep clean

White: Removes surface stains

