

PHILIPS

sonicare



Sonic electric
toothbrush with app

DiamondClean
Smart

Personalized coaching

Smart brush head sensors

Smart brush head recognition

5 modes, 3 intensities



HX9985/01

Our best ever toothbrush, for complete oral care

Feel confident, with 100% coverage

The DiamondClean Smart is our best toothbrush, for complete care^{*****}. High-performance brush heads let patients focus on all areas of their oral health, and our Smart Sensor technology gives personalized feedback and coaching.

Superior results for patients' oral care needs

- Removes 10x more plaque* for a deep clean
- Up to 7x healthier gums* in just 2 weeks
- Remove up to 5x more stains* in just 3 days
- Supports long-lasting clean breath

Personalized feedback with smart sensor technology

- Tracks and improves your patients' coverage
- Alerts patients when they're brushing too hard
- Guides patients to reduce their scrubbing
- A second chance to clean missed spots
- Helps patients start to take better care of trouble areas

Intuitive brushing experience

- Brush heads sync with optimal mode for ultimate results
- Offers 5 modes, 3 intensity settings
- Make sure patients get the most from their brush heads

Highlights

Say goodbye to plaque



C3 Premium Plaque Control is our deepest cleaning brush head yet. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact and remove up to 10x more plaque from hard-to-reach areas*.

Care for your gums



Click on the G3 Premium Gum Care brush head to improve your gum health. Its smaller size and targeted gum line bristles deliver gentle yet effective cleaning along the gum line, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation* and up to 7x healthier gums in just two weeks.*

W3 Premium White brush head




W3 Premium White is our best brush head for removing stains. With its densely packed central stain removal bristles, it's clinically proven to remove up to 5x more stains* in just three days.

TongueCare+ tongue brush



With 240 elastomer MicroBristles, the TongueCare+ tongue brush is specifically designed to clean the soft and porous surface of the tongue. The flexible MicroBristles contour to and gently clean around the ridges and grooves, to remove bacteria buildup and drive our bacteria-killing BreathRx tongue spray deeper. The compact shape allows patients to comfortably clean the entire tongue.


Focus areas and goal-setting



If you identify specific spots in a patient's mouth that need more attention due to plaque buildup, gum recession or any other issue, DiamondClean


Smart can help them tackle those trouble areas. Patients can highlight these spots in their personalized 3D mouth-map within the Philips Sonicare app, and they will receive reminders to give these areas the extra care they need each time they brush. The Philips Sonicare app maintains a running history of their brushing data, so patients can easily review their performance, track their progress and work to improve their daily oral care routine.

Location sensor



DiamondClean Smart's location sensor shows patients where they're brushing too little for better coverage. If there are spots your patients consistently miss when they're brushing, the location sensor will bring them to their attention.

Scrubbing sensor



With an electric toothbrush, patients should let the brush do the work rather than using a scrubbing motion. DiamondClean Smart has a scrubbing sensor that helps guide them to reduce scrubbing for optimized technique and a gentler clean.


DiamondClean Smart

Sonic electric toothbrush with app

HX9985/01


Highlights

BrushSync mode pairing



Whatever their oral care needs, DiamondClean Smart can help your patients achieve them. A microchip in their brush head tells their DiamondClean Smart which one they're using. So if they click on a gum-care brush head, for example, their toothbrush knows to select the optimal mode and intensity to provide a gentle yet effective care for their gums. All they need to do is to press the power button.

Smart replacement reminders



Over months of normal use, brush heads gradually wear down and lose their stiffness. An algorithm combining the data from smart brushing sensors and the microchip in the brush heads gives patients an accurate replacement reminder for each brush head. The Philips Sonicare app includes a smart brush head performance monitor to track brush head effectiveness based on patients' actual brushing time and pressure. It tells them when a brush head needs to be replaced so they continue to


get the best results. They even have the option to order replacement brush heads from within the app.

5 modes, 3 intensity settings




With five brushing modes (Clean, White+, Deep Clean+, Gum Health, TongueCare+) and three intensity settings (Low, Medium, High), your patients can personalize their brushing experience for optimal comfort and a complete clean. See Specifications for more detailed information on modes.

Pressure sensor



Your patients likely won't notice if they're brushing too hard, but their DiamondClean Smart will. If they apply too much pressure, the intuitive pressure sensor will make the toothbrush handle gently pulsate so they know exactly when to adopt a softer touch. In studies, 7 out of 10 aggressive brushers successfully reduced their brushing pressure thanks to the pressure sensor.

TouchUp feature



If your patients happen to miss any spots during a brushing session, DiamondClean Smart's TouchUp feature gives them the opportunity to go back for a second pass and immediately address missed areas to achieve a more complete clean.

Travel case and charging glass



It's easy for patients to take their DiamondClean Smart anywhere and everywhere with the new premium USB charging case. After placing their toothbrush in the case, they connect it to a laptop via USB or plug into a wall socket. Thanks to the charger code component, their charging cable will always be with them. Using the charger glass, patients simply place their Philips Sonicare electric toothbrush in the glass to charge. It can also be used to rinse after brushing. From a full charge, patients can enjoy two weeks of regular use.

DiamondClean Smart

Sonic electric toothbrush with app

HX9985/01

Specifications

Connectivity

Bluetooth® wireless technology: Connected brushing app

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): 14 days*****

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: White

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Display: Illuminated display

Battery indicator: Illuminated icon indicates battery life

Brush head recognition: Syncs with optimal mode

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Handle: 1 DiamondClean Smart

Brush heads: 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 W3 Premium White, 1 TongueCare+ tongue brush

Glass charger: 1

Travel case: USB charging travel case

Cleaning performance

Health benefits: Up to 7x healthier gums****

Performance: Removes up to 10x more plaque*

Whitening benefits: Up to 100% less stains*****

Speed: Up to 62000 brush movement/min

Timer: BrushPacer and SmarTimer

Pressure feedback: Ring lights up in purple, Vibrates handle to alert user

TouchUp: Ensures 100% coverage

Ingredients

Whitening toothpaste: 0.243% Sodium Fluoride, Zytex

Modes

3 intensities: High, Medium, Low

Clean: For exceptional everyday clean

Deep Clean+: For an invigorating deep clean

Gum Health: Special attention to molars

TongueCare+: For long lasting clean breath

White+: To remove surface stains

BrushSync Mode Pairing

C3 Premium Plaque Control: Pairs with Clean mode

G3 Premium Gum Care brush head: Pairs with Gum Health mode

W3 Premium White brush head: Pairs with White+ mode

TongueCare+ tongue brush: Pairs with TongueCare mode

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

Location sensor: Tracks & improves coverage

Scrubbing sensor: Guides to reduce scrubbing

3D mouth map: Highlights trouble areas

TouchUp: To address missed spots

