

Philips Sonicare
DiamondClean Smart 9300
Sonic electric toothbrush
with app

Personalized coaching

Smart brush head sensors
Smart brush head recognition
4 modes, 3 intensities



HX9903/21

Our best ever toothbrush, for complete oral care

Feel confident, with 100% coverage

The DiamondClean Smart is our best toothbrush, for complete care. Four high-performance brush heads let you focus on all areas of your oral health, and our Smart Sensor technology gives you personalized feedback and coaching.

Proven to improve oral health

- Removes up to 10x more plaque* for a deep clean
- Up to 7x healthier gums* in just 2 weeks
- Remove up to 100% more stains** in just 3 days

Personalized in-app coaching to perfect brushing

- Guides you to brush and highlights missed areas
- Takes you back to any missed spots for 100% coverage
- Start taking better care of your trouble areas
- Alerts you if you're pressing too hard
- Guides you to reduce your scrubbing

Optimize your brushing

- Make sure you're getting the most from your brush head
- Premium charging glass and travel case
- Brush heads automatically select the optimal settings
- Choose from 4 modes and 3 intensity settings

PHILIPS

sonicare

Highlights

Say goodbye to plaque



Click on the C3 Premium Plaque Control brush head to experience our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact*** and remove up to 10x more plaque from hard-to-reach areas*.

Care for your gums



Click on the G3 Premium Gum Care brush head to improve your gum health. Its smaller size and targeted gum line bristles deliver gentle yet effective cleaning along the gum line, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation* and up to 7x healthier gums in just two weeks.*


Whiter teeth, fast



Click on the W3 Premium White brush head to remove surface stains and reveal your whitest


smile. With its densely packed central stain removal bristles, it's clinically proven to remove up to 100% more stains in just three days.**

Never miss a spot




The Philips Sonicare app lets you know when you've achieved a thorough clean, and coaches you to be a better, more mindful brusher.

TouchUp feature



If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.


Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth


map, and remind you to pay extra attention to these areas.

Too much pressure?



You might not notice when you brush too hard, but your DiamondClean Smart will. If you use too much pressure, the light ring on the end of your handle will flash. This is a gentle reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

Easy does it



With an electric toothbrush, you let the brush do the work. We've built a scrubbing sensor into your handle as a gentle reminder to reduce scrubbing. This way, you can improve your technique and get a gentler, more effective clean.

Sonic electric toothbrush with app

Personalized coaching Smart brush head sensors, Smart brush head recognition, 4 modes, 3 intensities

HX9903/21

Specifications

Connectivity

- Bluetooth® wireless technology: Connected brushing app

Power

- Voltage: 110-220 V

Technical specifications

- Operating time (full to empty): 14 days*****
- Battery: Rechargeable
- Battery type: Lithium ION

Design and finishing

- Color: Pink

Service

- Warranty: 2-year limited warranty

Compatibility

- Android compatibility: Android phones, Bluetooth 4.0 enabled tablets
- iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher

Ease of use

- Battery indicator: Illuminated icon indicates battery life
- Display: Illuminated display
- Handle compatibility: Easy click-on brush heads
- Handle: Slim ergonomic design
- Brush head recognition: Syncs with optimal mode
- Replacement reminder: To always ensure best results, reminder icon lights up

Items included

- Travel case: Premium travel case

- Handle: 1 DiamondClean Smart
- Brush heads: 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 W3 Premium White
- Glass charger: 1

Cleaning performance

- Performance: Removes up to 10x more plaque*
- Health benefits: Up to 7x healthier gums****
- Speed: Up to 62000 brush movement/min
- Timer: BrushPacer and SmarTimer
- Whitening benefits: Up to 100% less stains*****
- Pressure feedback: Ring lights up in purple, Vibrates handle to alert user
- TouchUp: Ensures 100% coverage

Modes

- Clean: For exceptional everyday clean
- 3 intensities: High, Low, Medium
- Deep Clean+: For an invigorating deep clean
- Gum Health: Special attention to molars
- White+: To remove surface stains

BrushSync Mode Pairing

- C3 Premium Plaque Control: Pairs with Clean mode
- G3 Premium Gum Care brush head: Pairs with Gum Health mode
- W3 Premium White brush head: Pairs with White+ mode

Smart sensor technology

- Pressure sensor: Alerts when brushing too hard
- Location sensor: Tracks & improves coverage
- Scrubbing sensor: Guides to reduce scrubbing
- 3D mouth map: Highlights trouble areas
- TouchUp: To address missed spots
- Position sensor: Guides to areas you missed



Issue date 2022-01-26

Version: 2.2.3

EAN: 00 07502 00643 94

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Removes up to 7x more plaque vs. a manual toothbrush
** as compared to a manual toothbrush with a leading whitening toothpaste
** compared to DiamondClean
*** in Gum Health mode vs a manual toothbrush at 2 weeks
**** after 3 days of use in White+ mode vs a manual toothbrush
***** based on two periods of two-minute brushings per day on standard mode