



Philips Sonicare
HealthyWhite
Sonic electric toothbrush

3 modes
2 brush heads

HX6782/02

Whiter, healthier teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. HealthyWhite is proven to remove everyday stains in just two weeks, with the regular use of the Clean & White Mode.

Proven to improve oral health

- Clinically proven safe and gentle
- Naturally whiter teeth with patented sonic technology

Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

Offers a customized brushing experience

- Easy-start program to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2 minute mode
- Sensitive Mode: Gentle teeth and gum cleaning

PHILIPS
sonicare

Highlights

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Naturally whiter teeth



The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has proven to remove everyday stains, to help keep teeth naturally whiter.

Safe and gentle



Philips Sonicare is a gentle electric toothbrush for braces (brush heads wear out sooner when used on braces), and is safe for dental restorations (fillings, crowns, veneers) and periodontal pockets, too.

SmarTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Quadpacer




30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Clean and White Mode




2 minutes of Clean mode with an additional 30 seconds of White mode to focus on the visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Clean Mode

For exceptional everyday clean. Helps maintain whitening treatments.

Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums



Specifications

Design and finishing

- Color: Frost White

Power

- Voltage: 110-220 V

Technical specifications

- Operating time (full to empty): Up to 3 weeks**
- Battery: Rechargeable
- Battery type: Lithium ION

Service

- Warranty: 2-year limited warranty

Ease of use

- Handle: Slim ergonomic design, Rubber grip for easy handling
- Brush head system: Easy click-on brush heads
- Battery indicator: Light shows battery status
- Brushing time: Up to 3 weeks**

Items included

- Travel case: 1
- Charger: 1
- Handles: 1 HealthyWhite
- Brush heads: 1 ProResults standard, 1 ProResults compact

Cleaning performance

- Performance: Removes up to 90% of stains*
- Health benefits: Helps improve gum health
- Speed: Up to 62000 brush movement/min
- Whitening benefits: Whitens teeth up to 2 shades
- Timer: Quadpacer and SmarTimer

Modes

- Clean: For exceptional everyday clean
- Sensitive: Gentle teeth and gum cleaning
- Clean & White: 2 min Clean + 30 sec White



Issue date 2022-06-05

Version: 3.3.2

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* **based on two periods of two-minute brushings per day, on clean mode

* Removes up to 7x more plaque vs. a manual toothbrush