

PHILIPS

Replacement tip

VisaCare

Normal



SC6891/02



Replacement Tip Normal

For use only with Philips VisaCare microdermabrasion: used in the comfort of your home, the combined exfoliation and vacuum massage gently stimulate your skin's natural renewal processes and address the first signs of aging.

Microdermabrasion tip Normal

- Replace every 6 months
- Microdermabrasion tip for more intense treatment

Easy to replace exfoliating tip

- Easy click-on tip



asimpleswitch.com

Highlights

Microdermabrasion tip

Recommended to use for women with normal skin. Developed with dermatologists, safe for use at home. When using a new tip, always start again with a limited number of passes; allowing your skin to get used to iets more intense treatment.*

Easy fit

Easy to fit and compatible with Philips SC6220 and SC6240.

Exfoliation tip Normal

For optimal results, it is recommended to replace the microdermabrasion tip every six months.

Specifications

Benefits

Type of skin: For all skin types

Ease of use

Cleaning: Wash under fresh water, with soap and brush

Compatible with: Philips SC6220, Philips SC6240

Recommended replacement: Every 6 months

Replacement: Easy click-on tip



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

* compared to a tip used twice a week for six months.