



Philips RéAura  
Laser skin rejuvenation

SC5000/00

## Breakthrough anti-ageing technology

### Visible results

Philips RéAura laser skin rejuvenation is a breakthrough anti-ageing programme giving you radiant skin in 8 weeks. Based on Fraxel® laser technology, it stimulates your skin's natural cell renewal process, targeting several signs of ageing.

#### Exclusive advice from our Philips skincare experts

- Exclusive advice from the Philips skincare expert team

#### MyRéAura online treatment guide for best results

- Monitor your progress with your treatment diary on MyRéAura
- Personalise your treatment plan
- See what others have to say on the MyRéAura Community

#### Use Philips RéAura in three simple steps

- Step 3: Apply laser aftercare complex
- Step 1: Apply the laser performance gel
- Step 2: Apply the laser skin rejuvenation device

#### Visible results, confirmed by dermatologists\*

- Evens skin tone — 83%\*\*\*
- Reduces fine lines — 81%\*\*
- Smooths skin texture — 86%\*\*\*

# PHILIPS

# Highlights

## Reduces fine lines — 81%\*\*



Philips RéAura stimulates the production of collagen to reduce fine lines, for example around the eyes and mouth, and on the cheeks.\*\* Observed by dermatologists in 81% of users, after 12 weeks of home-use by 64 users. Results are based on analysis of before and after pictures by independent dermatologists.

## Evens skin tone — 83%\*\*



Philips RéAura achieves a more even skin tone and visibly diminishes age spots and sun spots.\*\*\* Observed by dermatologists in 83% of users, after 12 weeks of home-use by 64 users. Results are based on analysis of before and after pictures by independent dermatologists.

## Smooths skin texture — 86%\*\*




Philips RéAura encourages cell renewal in the deeper layers of the skin for smoother skin

texture.\*\*\* Observed by dermatologists in 86% of users, after 12 weeks of home-use by 64 users. Results are based on analysis of before and after pictures by independent dermatologists.

## Philips skincare expert team

Our skincare expert team has numerous years of beauty experience to share and is available just for you. They can help you to set up your personalised treatment plan, help interpret your skincare results or address concerns related to RéAura. You can contact them via phone, email and MyRéAura.

## Personalised treatment plan



MyRéAura allows you to define your unique treatment goals via a skin assessment. You can do this on [www.philips.co.uk/reaura](http://www.philips.co.uk/reaura), or have this done in the shop where you can buy RéAura. Our Philips skincare expert team is happy to help with your skin assessment. Just give them a call.

## See what others have to say

The MyRéAura Community is an exclusive place for you to connect with other RéAura users, share your experiences and see how others have experienced their treatment.

## Monitor your progress

Your MyRéAura diary is a great way to monitor the progress of your treatment plan and to track your results.

## Laser performance gel

The laser performance gel is a vital part of Philips RéAura laser skin rejuvenation. It provides optimal skin contact and smooth gliding of the Philips RéAura laser skin rejuvenation device. The gel is specially developed for Philips by the makers of Fraxel® laser technology. Without using gel, the skin may be treated unevenly. Dermatologist-tested. Hypoallergenic. Fragrance free. Non-comedogenic (will not clog pores).

## Laser skin rejuvenation device

During treatment, thousands of microscopic laser beams pass down into the epidermis and the dermis. This triggers a natural process of cell renewal. Old and damaged skin cells are replaced with new, fresh skin cells, containing more collagen. Repeated treatments boost the collagen network, making it denser. The result is a smoother skin texture and the reduction of fine lines. At the same time, pigmented skin cells are pushed out of the skin. This gives your skin a more even tone.

## Laser aftercare complex

The laser aftercare complex is a vital part of Philips RéAura laser skin rejuvenation. It hydrates the skin during the cell renewal process, after using the Philips RéAura laser skin rejuvenation device. The aftercare complex contains calming chamomile as well as St. John's Wort and edelweiss extracts, which are known for their anti-inflammatory properties. Formulated for Philips by the makers of Fraxel® laser technology. Dermatologist-tested. Hypoallergenic. Fragrance free. Non-comedogenic (will not clog pores).

# Specifications

## Items included

- Charger: Multi voltage charger
- Instructions for use: Demo DVD, Full user manual
- Storage: Pouch

## Application areas

- Face and body: Face, Neck, Chest, Hands, Arms

## Safety and adjustable settings

- 3 intensity settings: 1 to practise, 2 and 3 for results
- Child lock: Detachable tip as child lock
- Eye protection: Impossible to shine into eyes
- Overtreat protection: Locks for 24 hrs when overtreating

## Power

- Voltage: 100-240 V
- Charging time: 60 minutes, First charge 3 hours
- Power system: Rechargeable battery
- Running time: 1 charge treats all body areas

## Technical specifications

- Laser class: 1M
- Laser type: Diode
- Output wavelength: 1435 nm

## Service

- Warranty: 2 years global warranty



Issue date 2024-01-14

Version: 7.1.1

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* After 12 weeks of home use by 64 users. Results are based on analysis of before and after pictures by independent dermatologists.