


Philips VisaCare  
Replacement tip

Sensitive

SC6890/02



## Replacement Tip Sensitive

For use only with Philips VisaCare microdermabrasion: used in the comfort of your home, the combined exfoliation and vacuum massage gently stimulate your skin's natural renewal processes and address the first signs of aging.

### **Easy to replace exfoliating tip**

- Easy click-on tip

### **Microdermabrasion tip Sensitive**

- Microdermabrasion tip for mild treatment
- Replace every 6 months

**PHILIPS**

Replacement tip  
Sensitive

SC6890/02

# Specifications

## Benefits

- Type of skin: For sensitive skin

## Ease of use

- Compatible with: Philips SC6220, Philips SC6240

- Cleaning: Wash under fresh water, with soap and brush
- Replacement: Easy click-on tip
- Recommended replacement: Every 6 months

# Highlights

## Microdermabrasion Sensitive

Recommended to use for women new to microdermabrasion treatments, or for women with sensitive skin or in periods when the skin is more sensitive due to seasonal or hormonal changes. Developed with dermatologists, safe for home use. When using a new tip, always start again with a limited number of passes; allowing your skin to get used to a more intense treatment.\*

## Easy fit

Easy to fit and compatible with Philips SC6220 and SC6240.

## Replace every 6 months

For optimal results, it is recommended to replace the microdermabrasion tip every six months.



Issue date 2022-06-25

Version: 6.1.1

EAN: 08 71010 36752 59

© 2022 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* compared to a tip used twice a week for six months.