



Philips Sonicare AirFloss
Interdental - Rechargeable

w/ 1 nozzle



HX8211/02

Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

Improves oral health

- Removes more plaque between teeth
- Helps improve gum health in just two weeks
- Helps prevent cavities between teeth

Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

Easy to use

- Air and micro-droplet technology
- Comfortable ergonomic handle
- Guidance tip ensures correct placement

Helps build a healthy habit

- An easy way to improve cleaning between teeth

PHILIPS
sonicare


Highlights

Air and micro-droplet tech



AirFloss utilizes unique air and micro-droplet technology to send a burst of pressurized air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

An easy way to start




Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

30 Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

Guidance tip




The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

Removes more plaque between teeth



Removes more plaque between teeth

Ergonomic handle

Designed to be easy to maneuver in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.

Specifications

Power

- Voltage: Multi-voltage charger

Technical specifications

- Operating time: Replace nozzle every 6 months
- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Design specifications

- Color(s): White with green accents

Service

- Warranty: 2-year limited warranty

Ease of use

- Battery indicator: Flashing indicates low battery

- Handle: Slim ergonomic design
- Nozzle attachment: Easily snaps on and off
- Battery Life: Up to 3 wks when fully charged

Items included

- AirFloss Handle: 1
- Charger: 1 multi-voltage charger
- AirFloss Nozzle: 1
- Nozzle holder: 1

Cleaning performance

- Health benefits: Improves gum health in only two weeks
- Nozzles: Best results change every 6mos
- Auto-burst feature: Cleans entire mouth in 30 sec
- Performance: Removes up to 5x more plaque*



Issue date 2023-01-14

Version: 5.5.2

EAN: 08 71010 36233 28

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* than manual toothbrushing alone