



Philips Lumea
IPL hair removal system

For men

TT3003/61

Prevent hair regrowth

enjoy smooth skin everyday

Philips Lumea Plus IPL hair removal system TT3003 uses gentle pulses of light to stimulate the hair root into a resting phase. Hair sheds naturally and regrowth is inhibited. Using it every two weeks enables you to enjoy smooth skin everyday.

Effective hair regrowth control

- Body hair regrowth control

Intense Pulsed Light (IPL) technology

- Suitable for home use

Application areas

- Safe for full body use

Suitable hair and skin types

- NOT suitable for red, grey, light-blond hair & skin type VI

PHILIPS

Specifications

Technical specifications

- Voltage: 100-240 V
- Battery type: Lithium Ion
- Power: 7.5 W
- Frequency: 50-60 Hz
- Light energy level: 2 - 6.5 J/cm² (depending on chosen setting)
- Light spectrum: >570 nm
- High performance lamp: Delivers > 250,000 flashes

Features

- 5 light energy settings: Adjustable to your personal skin type
- integrated UV-filter: Protects your skin from UV light
- Ready-to-flash light: Indicates correct application on the skin
- 4cm² application window: For an effective treatment of large areas

Accessories

- Storage pouch: included

- Cleaning cloth: included
- AC/DC Adaptor/Charger: included

Power

- Voltage: 100-240 V
- Charging time: 60 minutes
- Power system: Rechargeable battery
- Running time: Minimum 370 flashes

Items included

- Instructions for use: Full user manual, Quick Start Guide

Service

- Warranty: 2 years global warranty

Safety and adjustable settings

- Safety ring system: Prevent unintentional flashing

Logistic data

- Country of origin: Hungary

Highlights


IPL technology

Philips Lumea works with an innovative light-based technology called "Intense Pulsed Light" (IPL) derived from technology used in professional beauty salons for hair removal. Philips Lumea now brings this innovative new technology, suitable for home use, to be used by you wherever and whenever you want.

Suitable hair and skin types

Philips Lumea works effectively on (naturally) dark-blond, brown and black hair. Philips Lumea is NOT effective for (naturally) red, grey and light blond hair. Philips Lumea is NOT suitable for very dark skin (skin type VI)

Application body areas



Philips IPL for men has been developed with leading dermatologists for safe, effective and gentle application on arms, underarms, shoulders, stomach, pubic area (excluding the scrotum), legs, chest and back.

Body hair regrowth control

Philips Lumea for Men applies gentle pulses of light to the hair root. As a consequence the hair sheds naturally and hair regrowth is inhibited. Repeating this treatment every two weeks for the underarms, pubic area, chest, stomach, legs and back will keep your skin smooth every day. For a most effective treatment pre-shave the application areas with the included Bodygroom as long as you still experience regrowth

