



Philips VisaCare
Microdermabrasion

Even Skin Tone
DualAction
Air Lift & Exfoliation system
For all skin types

SC6250/02

Flawless, even looking skin

with Philips VisaCare Prestige

Fade dark spots for flawless, even skin with VisaCare Prestige's salon-proven microdermabrasion technology. Its unique DualAction Air Lift & Skin exfoliation system has been specifically designed to deliver skin tone improvement.

Easy to use

- Designed for every day use
- Charge overnight for 3 weeks of use
- Cordless usage

Flawless, even looking skin

- Air lift massage for an immediate healthy glow
- Gentle skin exfoliation with 150.000 crystal particles

Safe and gentle on your skin

- Developed with expert scientists and dermatologists
- For all skin types and colors, from lightest to darkest

Gentle on your skin

- Can be used on wet skin for a smoother glide

PHILIPS

Specifications

Ease of use

- Cordless use
- Battery indicator: illuminated icon

Power

- Voltage: 100-240 V
- Charging time: 8 hours
- Power system: Rechargeable battery

Items included

- Cleaning brush: To clean VisaCare tips

- Power adapter: 100 - 240 V adapter
- Storage pouch: For storage and travel
- Exclusive and convenient stand: charging and storing stand
- Brightening/daily tip: can be used everyday
- Normal tip: twice a week for firming feeli

Service

- Warranty: 2-year limited warranty

Highlights

Gentle skin exfoliation

150,000 crystal particles gently sweep away the top layer of the skin where pigment is accumulated to reveal smoother, fresher skin.

Air lift massage

Air lift massage increases the blood flow by 67% giving skin an immediate healthy glow.

For use every day

VisaCare Even Skin Tone is so gentle on your skin that it can be used every day, together with your other skin products. The full treatment takes just 2 minutes.

Cordless usage

VisaCare Even Skin Tone is cordless so it's comfortable and convenient to use.

Overnight charge

Charge your VisaCare Even Skin Tone in the charging stand included. You'll get 3 weeks of treatment from just one overnight charge.

Developed with experts

We developed VisaCare Even Skin Tone together with expert scientists and dermatologists.

For all skin types

VisaCare Even Skin Tone can be used on all skin types, from the lightest to the darkest.

Use on wet & dry skin



The device has been developed for best results on wet skin. Simply moisten your skin with some water before use. It can also be used on dry skin and, of course, it's up to you to decide which you prefer. Both types of usage will make your skin feel firmer, but you'll find that the tip glides over wet skin more easily for a more pleasant and gentle experience.

