

Philips VisaCare
Replacement tip

Normal

SC6891/02



Replacement Tip Normal

For use only with Philips VisaCare microdermabrasion: used in the comfort of your home, the combined exfoliation and vacuum massage gently stimulate your skin's natural renewal processes and address the first signs of aging.

Easy to replace exfoliating tip

- Easy click-on tip

Microdermabrasion tip Normal

- Replace every 6 months
- Microdermabrasion tip for more intense treatment

PHILIPS

Replacement tip
Normal

SC6891/02

Specifications

Benefits

- Type of skin: For all skin types

Ease of use

- Compatible with: Philips SC6220, Philips SC6240

- Cleaning: Wash under fresh water, with soap and brush
- Replacement: Easy click-on tip
- Recommended replacement: Every 6 months

Highlights

Microdermabrasion tip

Recommended to use for women with normal skin. Developed with dermatologists, safe for use at home. When using a new tip, always start again with a limited number of passes; allowing your skin to get used to a more intense treatment.*

Easy fit

Easy to fit and compatible with Philips SC6220 and SC6240.

Exfoliation tip Normal

For optimal results, it is recommended to replace the microdermabrasion tip every six months.



Issue date 2022-06-25

Version: 5.2.1

EAN: 08 71010 36753 34

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* compared to a tip used twice a week for six months.