

PHILIPS

sonicare

Standard sonic
toothbrush heads

e-Series

2 pack

Standard size

Screw-on

HX7002/20



Classic power. Classic clean.

A fresh start to better oral health

The Philips Sonicare E-Series is recommended for those with a focus on effective plaque removal every day.

A brush head that fits multiple handles

- Screw-on system design for Essence, Elite, Xtreme, Advance

Proven to improve oral health

- Part of a better oral healthcare routine

Designed to optimise performance

- Philips Sonicare's advanced sonic technology
- Our classic brush head
- Reminder bristles ensure your most effective clean

Highlights

Our classic brush head




This Philips Sonicare toothbrush head features a contoured profile to naturally fit the shape of your teeth for a thorough clean.

Maximised sonic motion



Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

Screw-on brush head system




The Philips Sonicare E-Series brush heads screws on and off your brush handle for a secure fit and easy maintenance and cleaning.

Better oral healthcare



Like all authentic Philips Sonicare branded brush heads, this brush head is safe on teeth and gums. Each brush head has been quality tested for exceptional performance and durability.

Reminder bristles



At first glance it may not be obvious, but brush heads lose stiffness and gradually wear down over months of normal use. Our blue reminder bristles fade to white and help you recognise when it's time for a replacement. For optimal results, replace your brush head every three months.

Specifications

Design and finishing

Colour: White

Design and finishing

Bristle stiffness feel: Medium

Reminder bristles: Blue bristle colour fades away

Size: Standard

Compatibility

Brush head system: Screw-on

Suitable for these models: Essence+

Items included

Brush heads: 2 e-Series standard

Quality and performance

Replacement: Every 3 months

Tested: for optimal usage

