

# PHILIPS

## sonicare

Compact Sonicare  
toothbrush head

ProResults

1-pack

HX6021



## Focussed Clean

### Sonicare electric toothbrush head

The smaller brush head size of the Philips Sonicare electric toothbrush heads helps you easily maneuver the brush to focus on hard-to-reach areas.

#### **A brush head that fits multiple handles**

- Snap-on brush head

#### **Easy to maneuver to focus on hard-to-reach areas**

- Small brush head size helps you focus in tight places

#### **Designed to optimize performance**

- Reminder bristles ensure you are getting an effective clean

#### **Thorough clean feeling in your mouth**

- Brush head design maximizes sonic motion
- Contoured bristles fit the natural shape of your teeth
- Angled, radial-trimmed bristles expands the area covered

## Highlights

### Snap-on brush head



Easy brush head placement and handle cleaning

### Angled radial-trimmed bristles



The curved surface area created by angled, radial-trimmed bristles helps maintain contact with the surfaces of the teeth and gums as the brush head pivots, covering a larger surface area without increasing the size of the brush head.

### Contoured bristles



These Philips Sonicare electric toothbrush bristles are trimmed in a contoured pattern to fit the natural shape of teeth. The peaks and valleys along the length of the brush head are


ergonomically designed to follow teeth topography and fit in between teeth so you're less likely to miss places when you brush.

### Maximize cleaning performance



Thoughtful design and precision engineering ensure that when this Philips Sonicare replacement brushes are combined with the sonic motion and flexible brushing modes, you get a thorough clean.

### Reminder bristles



Reminder bristles let's you know when to replace the brush head. After three months of normal use bristles exhibit fatigue, and brush heads are less effective. Replace your brush heads every 3 months.

### Small brush head

The contoured bristle design on a smaller surface cleans as thoroughly as other Sonicare brush heads while allowing you to reach problem areas your dentist may have pointed out. This design helps you reach and focus on tight places where plaque and stains can hide.

## Specifications

### Technical specifications

**Operating time:** For optimal results, a new brush head is recommended every three months., For optimal results, a new brush head is recommended every three months.

### Weight and dimensions

**Dimensions brush head packaging:** 18.3 h x 6.3 w x 4.8 d cm

**Weight brush head packaging:** 0.026, 0.026 kg

### Items Included

Hygienic travel cap: 1 pcs

### Ease of use

**Suitable for these models:** FlexCare+, FlexCare, HealthyWhite, FlexCare+, HydroClean, FlexCare, EasyClean, HealthyWhite, HydroClean, EasyClean

