

# PHILIPS



Power Toothbrush  
with SenseIQ

Sonicare 9900  
Prestige

Our most advanced toothbrush  
SenseIQ senses, adapts, cares  
All-in-One brush head  
AI-powered Sonicare app



HX9992/22

## Exceptional care personalized to you with SenseIQ

Sonicare 9900 Prestige with SenseIQ technology is our most advanced power toothbrush. It enhances our proven sonic technology with intelligence, sensing and adapting for holistic care so you always get it right.

### **Designed around you**

- Sleek and compact design

### **Always get it right with SenseIQ**

- Adapts to you in real time when you brush too hard
- Real-time feedback on your brushing
- SenseIQ technology senses, adapts and cares
- Senses your brushing style
- Cares intuitively for teeth and gums

### **Personalized clean without compromise**

- All-in-One brush head for complete care
- Personalized guidance from the Philips Sonicare app
- Sonic brushing: the masterclass in effective cleaning

### **Premium accessories**

- Easy charging with the minimalist stand and base
- Prestige travel case with USB for go-anywhere style

# Sonicare 9900 Prestige

Power Toothbrush with SenseIQ

HX9992/22


## Highlights

### SenseIQ technology



As you brush, SenseIQ technology senses pressure, motion, coverage and more up to 100 times per second. It then adapts the intensity if you push too hard. The intelligence of SenseIQ takes the guesswork out of brushing modes, so you can experience effortless care and better brushing over time.

### Premium All-in-One brush head



Clean without compromise with the A3 Premium All-in-One brush head. Our all-time best has angled bristles to help remove up to 20x more plaque even in hard-to-reach areas\*. Triangular tips work to remove up to 100% more stains in less than two days\*\*\*. And longer bristles clean deep for gums that are up to 15x healthier in just two weeks\*\*. All that without ever switching brush heads. Did you know: brush heads become less effective after 3 months of use? Our BrushSync™ feature can remind you when it's almost time for a replacement.

### Personalized guidance



The Philips Sonicare app is powered by artificial intelligence and works in perfect harmony with the toothbrush. Use it to brush with real-time guidance about pressure, motion, coverage, duration and frequency. View progress reports by day, week, month and year. Get personalized recommendations and actionable steps for better brushing. Meanwhile, auto-syncing keeps brushing data up-to-date in the app, even when you don't brush with the app by your side.

### Senses your brushing style



While you clean, Philips Sonicare 9900 Prestige uses sensors to detect the pressure you apply, cleaning motion and coverage – up to 100 times per second. The sensors also keep an eye on how long and how often you brush.

### Adapts to you in real time



Most of us apply too much pressure during brushing. That's why our toothbrush automatically adjusts intensity when you push too hard. And with the free Philips Sonicare app,


you can experience tailor-made guidance and personalized recommendations to take your brushing to the next level.

### Cares for your oral health




Philips Sonicare Diamond Clean Prestige 9900 is designed around you. Its state-of-the-art technology is so in tune with you and effortless to use that you barely know it is there. Luckily, the results say it all.

### Proven sonic brushing



Brush away plaque with sonic technology – our benchmark in effective and gentle cleaning. While you brush, up to 62,000 bristle movements and dynamic fluid action help clean even hard-to-reach areas, leaving your mouth feeling exceptionally fresh and clean.

### Real-time feedback



You might not notice when you brush too hard or start scrubbing, but your toothbrush will. The light ring on the end of your handle will gently remind you to ease off the pressure or scrubbing motion.

# Sonicare 9900 Prestige

Power Toothbrush with SenseIQ

HX9992/22

## Highlights

### Prestige travel case with USB



The elegant and compact travel case can be slipped into any bag for carefree convenience while away from home. The durable case is wrapped in soft-touch vegan-friendly leather. And for extra ease, it includes an integrated USB-

C port, and a USB charging cable, so you can easily recharge the toothbrush without removing it from the case.

### Easy charging

When it is time to charge the toothbrush, simply place it on the charging base. Perfect for admiring your toothbrush's crisp, modern design.

### Sleek and compact



The Sonicare 9900 Prestige experience is about big technology delivered with style and compactness. The toothbrush is 18 mm shorter than DiamondClean and features premium materials to make this luxurious device easy on the eye, and pleasant to the touch.

# Sonicare 9900 Prestige

Power Toothbrush with SenseIQ

HX9992/22

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Power supply:** 110-220V

### Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** 2+ weeks\*\*\*\*

**Battery type:** Lithium ION

### Consumer Trade Item

**Net weight:** 1.58 kg

**Gross weight:** 1.984 kg

**GTIN:** 8710103995388

**Country of origin:** CN

### Outer Carton

**GTIN:** 18710103995385

### Design and finishing

**Color:** Midnight Blue

### Service

**Warranty:** 2-year limited warranty

### Pallet Layer

**Height:** 21 cm

**Width:** 101.6 cm

**Length:** 121.9 cm

### Compatibility

**Android compatibility:** Android OS 8.0 or above

**iOS compatibility:** iPhones on iOS 13.0 or above

### Ease of use

**Handle compatibility:** Easy click-on brush heads

**Handle:** Sleek and compact design, Seamless touch control

**Battery indicator:** Illuminated icon indicator

**Travel case:** Elegant, compact charging case, Vegan-friendly leather

**BrushSync Replacement Reminder:** Reminder icon lights up, to always ensure best results

### Items included

**Handle:** 1 Sonicare 9900 Prestige

**Brush heads:** 1 A3 Premium All-in-One, 1

TongueCare+ tongue brush

**Prestige USB charging case:** 1

**Charger base:** 1

**USB cable:** 1

### Cleaning performance

**Performance:** Removes up to 20x more plaque\*

**Health benefits:** Up to 15x healthier gums\*\*

**Speed:** Up to 62000 brush movement/min

**Timer:** BrushPacer and SmarTimer

**Whitening benefits:** Up to 100% less stains\*\*\*

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**Scrubbing feedback:** Ring lights up in amber, Vibrates handle to alert user, Option to turn on/off in app

### SenseIQ technology

**Pressure sensor:** Alerts when brushing too hard

**Scrubbing sensor:** Guides to reduce scrubbing

**Position sensor:** Tracks brushing coverage, for each brushing session

**3D mouth map:** Highlights missing areas

**Adaptive intensity:** Adjusts intensity level, when too much pressure applied

**Artificial Intelligence:** Detects and syncs brushing

### Sonicare App with A.I.

**Real time guidance:** Personalized recommendations, Day/month/year progress report

**Customizable modes:** Clean, White+, Gum health, Deep Clean, Sensitive

