




Philips Sonicare AirFloss  
Interdental - Rechargeable

w/ 1 nozzle



HX8211/02

## Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard-to-reach areas.

### Improves oral health

- Helps improve gum health in just two weeks
- Removes more plaque between teeth
- Helps prevent cavities between teeth

### Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

### Easy to use

- Air and micro-droplet technology
- Comfortable ergonomic handle
- Guidance tip ensures correct placement

### Helps build a healthy habit

- An easy way to improve cleaning between teeth

**PHILIPS**  
**sonicare**


# Highlights

## Air and micro-droplet tech



AirFloss utilises unique air and micro-droplet technology to send a burst of pressurised air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

## An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

## 30-Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

## Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

## Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

## Guidance tip




The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

## Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

## Removes more plaque between teeth



Removes more plaque between teeth

## Ergonomic handle

Designed to be easy to manoeuvre in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.

# Specifications

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Operating time: Replace nozzle every 6 months
- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

## Design specifications

- Colour(s): White with green accents

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Handle: Slim ergonomic design

- Battery indicator: Flashing indicates low battery
- Nozzle attachment: Easily snaps on and off
- Battery Life: Up to 3 weeks when fully charged

## Items included

- AirFloss Handle: 1
- Charger: 1 multi-voltage charger
- AirFloss Nozzle: 1
- Nozzle holder: 1

## Cleaning performance

- Performance: Removes up to 5x more plaque\*
- Nozzles: Best results change every 6m
- Auto-burst feature: Cleans entire mouth in 30 secs
- Health benefits: Improves gum health in only two weeks



Issue date 2023-01-14

Version: 5.0.2

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* than manual tooth-brushing alone