


Philips Sonicare
AirFloss Pro/Ultra -
Interdental cleaner

1 nozzle



HX8431/02



Healthier gums in 2 weeks*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra**** is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.**

Improves oral health

- Clinically proven as effective as floss for gum health**
- Up to 99.9% plaque removal***
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings


The easiest way to effectively clean between teeth

- Takes just 60 seconds: Point, Press, Clean!
- Total daily confidence of a fresh and clean mouth
- An easy way to start a healthy habit

PHILIPS
sonicare

Highlights

Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.***

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Fresh & clean with mouthwash




Fill the reservoir on the handle with either mouthwash or water, then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

Triple burst customization




AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Specifications

Design and finishing

- Color: Pink

Power

- Voltage: Multi-voltage charger

Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Charge time: 24 hours to full charge

Service

- Warranty: 2-year limited warranty

Ease of use

- Battery indicator: Handle LEDs show charge level
- Handle: Slim ergonomic design

- Nozzle attachment: Easily snaps on and off
- Battery Life: 1-3 weeks****

Items included

- Charger: 1
- AirFloss Pro/Ultra handle: 1
- AirFloss Pro/Ultra nozzle: 1

Cleaning performance

- Health benefits: Cleans entire mouth in 60 sec, Improve gum health in 2 weeks*
- Performance: Removes up to 99.9% plaque**
- Clinically proven results: As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction
- Customizable settings: Auto-burst or manual burst, Single/double/triple burst***
- For best results: Change nozzle every 6 months



Issue date 2023-03-16

Version: 2.2.1

© 2023 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.
*** From the treated areas; In a lab study, actual in-mouth results may vary
**** AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.
***** depending on burst setting used