

Philips Sonicare AirFloss  
Interdental - Rechargeable

**Rechargeable**  
w/ 1 nozzle



HX8111/02

## Philips AirFloss - An easier way to floss

Gently and effectively improves oral health

If flossing isn't a regular part of your oral care routine, you're not reaching the bacteria that can build up between teeth. The Philips Sonicare AirFloss helps you to get a deep clean everyday by taking the hassle out of flossing.

### **Achieve a 30-second clean**

- Simple one-button use; as easy as a mouse-click

### **Proven improvement of oral health**

- Up to 99.9% plaque removal

### **Simple: point & press**

- Easy-to-fill reservoir
- Guidance tip for easy placement
- Microburst technology cleans at the touch of a button
- Slim, angled nozzle for easy to even the hard-to-reach areas
- Slim, easy to hold handle
- Fits standard Sonicare chargers
- Charging after two weeks

**PHILIPS**  
**sonicare**

# Highlights

## Microburst technology



Microburst technology applies a quick burst of pressurized air and micro-water droplets to clean deep between teeth where a toothbrush can't reach

## Guidance tip



Simply glide the guidance tip along the gum-line until you feel it settle between the teeth. Provides easy positioning even on those hard-to-reach back teeth.

## Fits Sonicare chargers



You only need one charger to charge your Sonicare products.

## Simple one-button use




Simple one-button use for easy cleaning in 60 seconds.

## Easy-to-fill reservoir



Easy-to-fill reservoir uses only a tiny amount of water for less mess. Fill with mouthwash for enhanced oral health benefits.

## Up to 99.9% Plaque Removal



AirFloss uses microdroplets of air and water or mouthwash to remove up to 99% more plaque in-between teeth\*

## Slim, angled nozzle



Slim, angled nozzle for easy to even the hard-to-reach areas

## Slim, easy to hold handle



Easy to manoeuvre around the mouth.

## Charging after two weeks

Philips AirFloss lasts two weeks between charging.

# Specifications

## Power

- Voltage: Multi-voltage charger

## Technical specifications

- Operating time: Replace nozzle every 6 months
- Battery: Rechargeable
- Operating time (full to empty): 2 weeks or 14 full-mouth uses
- Charging time: 24 hour(s)

## Service

- Warranty: 2-year limited warranty

## Ease of use

- Battery indicator: Light in power button blinks to indicate when to recharge
- Handle: Slim ergonomic design



Issue date 2023-03-05

Version: 2.2.1

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* Compared to manual toothbrushing alone