www.philips.com/One

©2021 Koninklijke Philips N.V. (KPNV). All rights reserved. Philips and the Philips shield are trademarks of KPNV. Sonicare and the Sonicare logo are trademarks of Philips Oral Healthcare, LLC and/or KPNV.

©2021 Koninklijke Philips N.V. (KPNV). 保留所有权利。飞利浦及飞利浦盾牌标志是 KPNV 的商标。Sonicare 及 Sonicare 徽标是 Philips Oral Healthcare, LLC 和/或 KPNV 的商标。

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English

### Important safety information

Only use the product for its intended purpose. Read this important information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries.

### Warnings

- Make sure the handle is completely dry before you connect the USB-C cable.
- This appliance can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision.
- Children shall not play with the appliance.
- Only use original Philips accessories or consumables. Only use the USB cable provided with the product.
- This appliance contains no user-serviceable parts. If the appliance is damaged, stop using it and contact the Consumer Care Center in your country (see 'Warranty and support').
- Do not charge the product outdoors or near heated surfaces.
- Do not clean any part of the product in the dishwasher.
- This appliance has only been designed for cleaning teeth, gums and tongue.
- Stop using the appliance and consult your dentist/doctor if excessive bleeding occurs after use, if bleeding continues to occur after 1 week of use or if you experience discomfort or pain.
- If you have had oral or gum surgery in the previous 2 months, consult your dentist before you use this appliance.
- If you have a pacemaker or other implanted device, contact your doctor or the manufacturer of the implanted device prior to use.
- If you have medical concerns, consult your doctor before you use this appliance.
- This appliance is a personal care device and is not intended for use on multiple patients in a dental practice or institution.
- Stop using a brush head with crushed or bent bristles. Replace the brush head every 3 months or sooner if signs of wear appear.
- If your toothpaste contains peroxide, baking soda or bicarbonate (common in whitening toothpastes), thoroughly clean the brush head with soap and water after each use. This prevents possible cracking of the plastic.
- Avoid direct contact with products that contain essential oils or coconut oil. Contact may result in bristles dislodging.
- Charge, use and store the product at a temperature between 0 °C and 40 °C.
- Keep product and batteries away from fire and do not expose them to direct sunlight or high temperatures.
- If the product becomes abnormally hot or smelly, changes color or if charging takes longer than usual, stop using and charging the product and contact Philips.
- Do not place products and their batteries in microwave ovens or on induction cookers.
- Do not open, modify, pierce, damage or dismantle the product or battery to prevent batteries from heating up or releasing toxic or hazardous substances. Do not short-circuit, overcharge or reverse charge batteries.

- If batteries are damaged or leaking, avoid contact with the skin or eyes. If this occurs, immediately rinse well with water and seek medical care.

### Electromagnetic fields (EMF)

This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

### Intended use

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at [www.philips.com/welcome](http://www.philips.com/welcome).

### Intended use

Philips One power toothbrushes are intended to remove adherent plaque and food debris from the teeth to reduce tooth decay and improve and maintain oral health. Philips One power toothbrushes are intended for consumer home use. Use by children should be with adult supervision.

### Your Philips One (Fig. 1)

- 1 Travel case lid
- 2 Travel case
- 3 Brush head
- 4 Handle
- 5 Power Button
- 6 Battery light indicator
- 7 End cap
- 8 USB-C charging port
- 9 USB-C charging cable

### Getting started

#### Attaching the brush head

- 1 Align the brush head so the bristles point in the same direction as the front of the handle. Press the brush head down onto the metal shaft (Fig. 2).

Note: It is normal to see a slight gap between the brush head and the handle (Fig. 3). This allows the brush head to vibrate properly.

#### Charging your Philips One

- 1 Open the end cap at the bottom of the toothbrush handle (Fig. 4) and plug the USB-C charging cable into the USB-C charging port (Fig. 5).
- 2 Connect the USB cable into a 5V USB wall adapter (Type A).

Note: Wall adapter not included. Only use 5V IPX-4 adapter (Type A).

- 3 Your Philips One beeps twice to indicate that the toothbrush is charging.
- 4 Your toothbrush continues to charge until it is fully charged.
  - When the toothbrush is charging, the battery indicator fades up and down in white.
  - When fully charged, the battery indicator lights up solid white and turns off after 30 seconds of being fully charged.

Note: It can take up to 8 hours to fully charge the battery, but you can use the toothbrush before it is fully charged.

When not charging your Philips One, be sure the end cap is securely closed.

#### Travel case

Your Philips One comes with a travel case to store your power toothbrush when traveling or between uses.

Put your Philips One toothbrush in the travel case with the bristles facing downwards. Snap close the lid of the travel case.

- To remove your Philips One toothbrush in the travel case, follow these steps:
- 1 Flip open the travel case lid.
  - 2 While holding both ends of the travel case, pinch the base of the toothbrush (Fig. 6).
  - 3 The toothbrush will raise. Lift the toothbrush.

#### Using your Philips One toothbrush

If you are using this toothbrush for the first time, it is normal to feel more vibration than when using a manual toothbrush. The Philips One uses microvibrations to help gently clean your teeth. Please follow the brushing instructions below for the best experience. Philips One comes with a SmarTimer and QuadPacer to help you brush evenly throughout your mouth.

**SmarTimer**  
The SmarTimer indicates that your 2-minute brushing cycle is complete by beeping 3 times and automatically turning off.

**Note:** Dental professionals recommend brushing no less than 2 minutes twice a day.

#### QuadPacer

The QuadPacer is an interval timer that beeps every 30 seconds to remind you to brush all teeth evenly by dividing your mouth into 4 sections (Fig. 7). Begin brushing in section 1 (outside upper teeth) and brush for 30 seconds before you move to section 2 (inside upper teeth). Continue brushing in section 3 (outside lower teeth) and brush for 30 seconds before you move to section 4 (inside lower teeth).

**Note:** If you press the power button after you started the brushing cycle, the toothbrush will pause. After a pause of 30 seconds, the SmarTimer and QuadPacer reset.

### Brushing instructions

Please follow the brushing instructions for the best experience.

- 1 Wet the bristles and apply a small amount of toothpaste.
- 2 Place the toothbrush bristles against the teeth at a slight angle (45 degrees) (Fig. 8), pressing gently to make the bristles reach the gumline or slightly beneath the gumline.

Note: The bristles should slightly flare.

- 3 Press the power button to turn on the toothbrush.
- 4 Brush your teeth with back and forth strokes so that the bristles reach all surfaces of the teeth. Change to another section when the QuadPacer beeps. Continue this motion throughout the 2-minute brushing cycle.
- 5 To clean the inside surfaces of the front teeth, tilt the brush vertically and make several up-and-down strokes (Fig. 9).
- 6 After you have completed the brushing cycle, you can spend additional time brushing the chewing surfaces of your teeth and areas where staining occurs (Fig. 10). You can also brush your tongue, with the toothbrush turned on or off, as you prefer.

Your power toothbrush should be safe to use on braces (brush heads wear out sooner when used on braces) and dental restorations (fillings, crowns, veneers) if they are properly adhered and not compromised. If a problem occurs, please follow up with a dental professional.

### Battery status

This Philips One toothbrush is designed to provide at least 60 brushing sessions after it has been fully charged for 8 hours. This allows sufficient charge for 1 month's usage if you brush 2 times a day (each session lasting 2 minutes). When the battery light indicator flashes amber it is time to charge your Philips One (see chapter 'Charging your Philips One'). The table below specifies the battery status:

Number of brushings remaining	Battery indicator	Beep sound
40-60 brushings	White- solid	-
7-39 brushings	White- fades up and down	-
≤6 brushings	Amber- flashing	5 beeps after brushing
0 brushings remain	Amber- flashing	5 beeps x 2

Note: When the battery is completely empty, the toothbrush turns off. The toothbrush will not turn on until battery is recharged. Toothbrush does not turn on while charging.

Note: Battery usage time is expected to reduce over product lifetime.

### Cleaning

You should clean your toothbrush and accessories regularly to remove toothpaste and other residue. Failure to clean your product may result in an unhygienic product and damage to it may occur.

**Warning: Do not clean product or accessories in dishwasher, microwave, with chemicals or with boiling hot water.**

#### Toothbrush handle

- 1 Remove the brush head and rinse the metal shaft area with warm water (Fig. 11). Make sure you remove any residual toothpaste.
- 2 Wipe the entire surface of the handle with a damp cloth.

Note: Do not tap the handle on the sink to remove excess water.

#### Cleaning brush head

- 1 Rinse the brush head and bristles after each use.
- 2 Remove the brush head from the handle and rinse the brush head connection with warm water at least once a week.

#### Cleaning travel case

Rinse the travel case weekly with warm water and wipe clean with a damp cloth.

#### Cleaning USB charging cable and port

- 1 Use compressed air to clean the USB port on the toothbrush handle and connector on the USB charging cable.
- 2 Use a slightly dampened cloth to wipe the cable.

Note: Avoid contact of the USB port and connections with water. If contact occurs, dry immediately.

### Storage

If you are not going to use the product for a long time, clean it (see 'Cleaning'), charge it fully, and store it in a cool and dry place away from direct sunlight.

### Replacement

#### Brush head

Replace the brush heads every 3 months to achieve optimal results. Only Philips One brush heads are compatible with Philips One.

### Locating the model number

Remove brush head from handle. Look at the top of the handle near the metal shaft for the model number (HY120X).

### Troubleshooting

This chapter summarizes the most common problems you could encounter with the appliance. If you are unable to solve the problem with the information below, visit [www.philips.com/support](http://www.philips.com/support) for a list of frequently asked questions or contact the Consumer Care Center in your country.

Problem	Possible cause	Solution
I cannot fully attach the brush head. There is a gap between the brush head and handle	This is normal. A small gap is necessary to allow movement in order to deliver the vibrations needed to clean your teeth	
My Philips One toothbrush vibration is less powerful than before	The battery is running low	Charge the battery (see 'Charging your Philips One').
The toothbrush stops working	The battery is empty	Charge the battery (see 'Charging your Philips One').

### Warranty and support

If you need information or support, please visit [www.philips.com/support](http://www.philips.com/support)

**Warranty**  
Philips warrants its products for a period of one year after the date of purchase. If any defect due to faulty materials and workmanship occurs within this one-year warranty period, Philips will repair or replace the product at its expense. Philips will only pay for repair or replacement if convincing proof is provided, for instance by means of a receipt, that the day on which service is claimed is within the warranty period. The warranty does not cover products and/or product parts that are subject to wear, that can be considered as consumable parts by their nature. The warranty is not valid if a defect is due to damage caused by incorrect use, poor maintenance (e.g. parts blocked by scale) or if alterations or repairs have been carried out by persons not authorized by Philips. For the appliance to be used correctly, the user should strictly adhere to all instructions included in the user manual and should abstain from any actions or uses that are described as undesired or which are warned against in the user manual.

#### Important warranty restrictions for certain products and materials

For some products the warranty may be subject to certain restrictions, for instance as a result of the type of material used. Warranty restrictions, if any, can be found under the heading 'Warranty restrictions' in chapter 'Warranty and support' of the user manual.

#### Support

Support during and after the warranty period can be obtained in all countries where the product is officially distributed by Philips. In countries where the product is not distributed by Philips, service is provided by the national Philips organization. In this case there may be some delay if the required spare parts are not readily available. Turn to your dealer or an authorized Philips service center for service. Contact the Consumer Care Center in your country if you need more information or if you have problems obtaining service. If there is no Consumer Care Center in your country, you can contact the Service Department of Philips.  
**China**  
Philips Consumer Care centre  
Building A1, #718 Lingshi Road,  
JingAn District, Shanghai  
4008 800 008  
**Taiwan**  
0800 231 099  
[www.philips.com/support](http://www.philips.com/support)

#### Warranty restrictions

The terms of the international warranty do not cover the following:  
- Brush heads.  
- Damage caused by use of unauthorized replacement parts.  
- Damage caused by misuse, abuse, neglect, alterations or unauthorized repair.  
- Normal wear and tear, including chips, scratches, abrasions, discoloration or fading.

### Recycling

This symbol means that electrical products and batteries shall not be disposed of with normal household waste.  
- Follow your country's rules for the separate collection of electrical products and batteries.

**Removal of built-in rechargeable battery**  
The built-in rechargeable battery must be

removed by a qualified professional when the product is discarded. Instructions for removal of built-in rechargeable batteries can be found on [www.philips.com/support](http://www.philips.com/support)

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